D-dimer

Instructions For Use REF 5250

**D-dimer
Fiche technique**

**D-dimer
Anleitung**

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Contents

English	1
Français	6
Deutsch	12
Italiano	17
Español	22

Latex Agglutination Test for Fibrin D-dimer.**Reagents for 80 determinations****For in vitro diagnostic use only.****I. INTENDED USE**

Helena D-dimer is a latex agglutination test for semi-quantitative determination of fibrin D-dimer.

II. PRINCIPLE

D-dimer containing moieties are formed by plasmin degradation of Factor XIIIa cross-linked fibrin. Elevated levels of D-dimer are found in clinical conditions such as deep vein thrombosis (DVT), pulmonary embolism (PE) and disseminated intravascular coagulation (DIC)^{1,3}. D-dimer levels rise during pregnancy and high levels are associated with complications². Helena D-dimer utilizes a monoclonal antibody reacting with fibrin D-dimer or fragment D of fibrin but not with intact fibrinogen⁶, allowing D-dimer determination in human plasma. Serum samples suited for FDP analysis can also be used.

No international standard for D-dimer exists at present. For numerical values Helena has aligned the Helena D-dimer with commercially available enzyme-linked immunosorbent assay (ELISA) tests.

III. REAGENTS**A. Reagent description****1. Helena D-dimer Latex**

1 vial containing 1.7 ml suspension of latex beads. Beads are coated with anti D-dimer monoclonal antibody.

2. Saline Solution

2 vials containing 8 ml each of buffered saline, pH 7.3. Contains 0.2 g/L sodium azide.

3. Positive Control Plasma

1 vial containing 1.0 ml lyophilized human plasma enriched with fibrin D-dimer.

4. Negative Control Plasma

1 vial containing 1.0 ml lyophilized human plasma.

B. Reagent preparation

1. Allow all vials to warm to room temperature before use, minimum 10 minutes.
2. Add 1.0 ml (1000 μ l) Saline Solution to each of the negative and positive control plasma vials. Control plasma may be stored at 2...6°C for one month. See XIII. Notes point 3.
3. Agitate the latex by repeatedly inverting the vial for 5 seconds immediately before use. If latex mass collects at the neck of the vial or is visible in the latex suspension, vortex vial for 10 seconds, to ensure any settled material is resuspended.

IV. WARNINGS AND PRECAUTIONS

POTENTIAL BIOHAZARD MATERIAL. Negative and Positive Control Plasma are of human origin. Each donor unit of source plasma used in these products has been tested and found negative for Hepatitis B antigens, HIV I and II antibodies, Hepatitis C antibodies, syphilis antibodies and H.T.L.V. I/II antibodies by FDA approved methods. Because no test method can offer complete assurance that no infectious agents are present, these plasmas should be treated at the Biosafety Level 2 as recommended for any potentially infectious human serum or blood specimen in the Centers for Disease Control/National Institutes of Health manual, "Biosafety in Microbiological and Biomedical

Laboratories", 1984. All wastes containing biological material should be properly labelled and stored separately from other wastes. Dispose of all waste materials according to prescribed international, national and local regulations.

SODIUM AZIDE. Both the Saline Solution and the D-dimer Latex contain sodium azide, which may react with lead and copper plumbing to form highly explosive metal azides. Materials discarded into a sink should be flushed with a large volume of water to prevent azide build-up.

The test should be used in conjunction with clinical observations and results of other laboratory tests.

V. STORAGE AND STABILITY

The reagents should be stored at 2...8°C and used prior to the expiration date recorded on the box. Control plasmas, once reconstituted may be stored at 2...6°C for one month.

VI. SAMPLE COLLECTION

Plasma anticoagulated with citrate is recommended. Collect samples in nine volumes of blood and one volume of 3.2% buffered sodium citrate (0.105 M). Centrifuge at 3000 x g for 10 minutes. Citrated plasma samples may be stored at room temperature for 2 hours or at 2...8° C for 18 hours. A single freeze-thaw cycle does not affect the assay response. Refer to NCCLS H21-A2 for specific sample collection guidelines.

VII. PROCEDURE

A. Materials provided

- 1 x Helena D-dimer Latex
- 2 x Saline Solution
- 1 x Positive Control Plasma
- 1 x Negative Control Plasma
- 50 Mixing sticks
- 16 Test cards for 6 samples each.

B. Materials required but not provided

- Pipettes for 20, 100 and 1000 μ l volumes
- Pipette tips
- 3 ml test tubes
- Timer

C. Qualitative method

1. Place 20 μ l of sample, positive and negative control plasmas in circles on a test card.
2. Place 20 μ l of latex suspension in a nearby area of each circle.
3. Quickly mix the sample and latex using clean mixing sticks for each sample. Start the timer.
4. Rock the test card gently and read agglutination between 180 and 200 seconds. Positive (+) or negative (-) agglutination are compared to results obtained using the controls. The positive control is merely qualitative and should not be further diluted. Non-agglutinated latex means the sample is normal and no further testing is required.

D. Semi-quantitative method

Performed only on samples tested positive:

1. Serially dilute 100 μ l of sample 1:2, 1:4 and 1:8 with 100 μ l Saline Solution using small test tubes.
2. Mark the positions of sample dilutions on the test card and mix with latex suspension as per VII.
 - C. The D-dimer concentration may be determined from the table in IX. EXPECTED VALUES.

VIII. QUALITY CONTROL

The positive and negative controls provided should be used for quality control of each kit. It is recommended that both positive and negative controls are tested each time the kit is used. If either positive or negative control fails to elicit the appropriate response, patient results obtained on that occasion should not be used and the kit should be discarded or returned to Helena.

IX. EXPECTED VALUES

Agglutination occurs within 180 to 200 seconds for samples containing more than 250 ng/ml D-dimer. By testing serially diluted plasmas, semi-quantitative results can be obtained.

D-dimer Level (ng/ml)

(ng/ml)	Undiluted	1:2	1:4	1:8
< 250	-	-	-	-
250-500	+	-	-	-
500-1000	+	+	-	-
1000-2000	+	+	+	-
> 2000	+	+	+	+

Agglutination may be more pronounced and appears more rapidly at higher D-dimer concentrations. If you wish to express the results in fibrinogen equivalent units (FEU) multiply the D-dimer levels in the table above by 2, e. g. <250 ng/ml becomes <500 FEU. As discussed in several reports⁴, some commercial latex test do not have the claimed sensitivity when compared to commercial ELISA.

X. INTERPRETATION OF RESULTS

Agglutination occurs within 180-200 seconds for samples containing more than 250 ng/ml D-dimer. The mean level of D-dimer in a healthy population is between about 8 and 135 ng/ml, and neat plasma from normal, healthy individuals should not agglutinate. The negative predictive value of the Helena D-dimer for thrombosis is high⁵. The circulatory half-life of D-dimer is about 12 hours. Elevated D-dimer levels can therefore persist for some time after the active process has ceased.

In clinical studies on normal subjects, patients with phlebographically confirmed DVT, patients with DIC and patients with pre-eclampsia (Pre-EC) the following results were obtained:

Helena D-dimer results†

Samples	η	-	1:1	1:2	1:4	≥1:8
Normal	101	100	1*	-	-	-
DVT	48	3	10*	7*	14*	14*
DIC	29	0	3	3	4	19
Pre-EC	6	2	1	3	-	-

† The symbols η and - denote the number of patients and a negative Helena D-dimer result, respectively. Titters indicate the highest dilution at which the sample shows agglutination.

* The agglutination was inhibited by addition of the D-dimer specific antibody (0.2 mg/ml) MA-8D3 but not with a non-related antibody PAM-1.

XI. PERFORMANCE CHARACTERISTICS

A. Specificity

The monoclonal antibody used in this device, MA 8D3, is specific for D-dimer by virtue of the screening method used for hybridoma selection³. A hybridoma, secreting antibodies which reacted positively with purified D-dimer but not with whole fibrinogen or fragment D of fibrinogen, was selected. No cross-reactivity with fibrinogen or des-AA-fibrinogen was observed when either analyte was substituted for plasma in this assay. Plasma from 16 patients with rheumatoid arthritis were assayed and 14 of these were found to be non-agglutinating with the Helena D-dimer test. The two agglutinations could be inhibited by addition of the D dimer specific monoclonal antibody MA-8D3 but not with addition of a monoclonal of the same subgroup, IgG1k (PAM 1). This suggests that the Helena D-dimer is insensitive to rheumatoid factor disturbances.

B. Reproducibility

Three plasma samples were selected to test reproducibility of the assay. Each sample was tested 10 times on each of three different days. In each case, where as positive result was obtained, the sample was titrated. The results in D-dimer (ng/ml):

Sample	D-dimer level	Result
Normal	<250 ng/ml	negative
Intermediate	3000 ng/ml	titer 1:8
High	> 16,000 ng/ml	titer 1:64

C. Accuracy

The Helena D-dimer kit was compared to another commercially available D-dimer latex. Both products gave a negative reaction when tested on 25 normal specimens. When 30 patient plasma samples were tested with ELISA and Helena D-dimer, all samples with ELISA values above 225ng/ml showed agglutination with Helena D-dimer.

XII.LIMITATIONS OF TEST

1. A negative D-dimer test does not completely rule out thrombosis. A negative predictive value for patients with suspected DVT has been found to be 94% with the Helena D-dimer kit⁵. Detection of elevated levels of D-dimer should be used with other clinical information in forming a diagnosis.
2. Agglutination with samples containing normal D-dimer levels may occasionally occur due to non-specificity.
3. A small number of samples, when mixed with the latex, may exhibit white flakes which should not be confused with agglutination.

XIII.NOTES

1. The results may be reported either in D-dimer units or in fibrinogen equivalent units (FEU); 1 ng/ml of D-dimer is equivalent to approximately 2 ng/ml of FEU.
2. Agglutination will be more pronounced and appears more rapidly at higher D-dimer concentrations.
3. The rubber stopper in the control vials may be discarded after vial opening.

XIV.REFERENCES

1. Elms, M.J. et al.: Rapid detection of cross-linked fibrin degradation products in plasma using monoclonal antibody-coated latex particles. *J. Clin. Path.* 35, 360-364, 1986.
2. Ballegeer, V. et al.: Fibrinolytic response to venous occlusion and fibrin fragment D-dimer levels in normal and complicated pregnancy. *Thromb. Haemostas.* 58, 1030-1032, 1987.
3. Declerck, P.V. et al.: Fibrinolytic response and fibrin fragment D dimer in patients with deep vein thrombosis. *Thromb. Haemostas.* 58, 1025-1029, 1987.
4. Bounameaux, H. et al.: Measurement of plasma D-dimer for diagnosis of deep venous thrombosis. *Am. J Clin. Path.* 91, 82-85, 1989.
5. Hansson P.O. et al.: Can laboratory testing improve screening strategies for deep vein thrombosis at an emergency unit? *J. Intern. Med.* 235, 143-151, 1994.
6. Holovet P et al.: Binding properties of monoclonal antibodies against human fragment D-dimer of cross-linked fibrin to human plasma clots in an in vivo model in rabbits. *Thrombosis and Haemostasis* 61(2), 307-313, 1989.

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1. Elms, M.J. et al.: Rapid detection of cross-linked fibrin degradation products in plasma using monoclonal antibody-coated latex particles. J. Clin. Path. 85, 360-364, 1986.
2. Ballegeer, V. et al.: Fibrinolytic response to venous occlusion and fibrin fragment D-dimer levels in normal and complicated pregnancy. Thromb. Haemostas. 58, 1030-1032, 1987.
3. Declerck, P.V. et al.: Fibrinolytic response and fibrin fragment D-dimer in patients with deep vein thrombosis. Thromb. Haemostas. 58, 1025-1029, 1987.
4. Bounameaux, H. et al.: Measurement of plasma D-dimer for diagnosis of deep venous thrombosis. Am. J Clin. Path. 91, 82-85, 1989.
5. Hansson P.O. et al.: Can laboratory testing improve screening strategies for deep vein thrombosis at an emergency unit? J. Intern. Med. 235, 143-151, 1994.

6. Holovet P et al.: Binding properties of monoclonal antibodies against human fragment D-dimer of cross-linked fibrin to human plasma clots in an in vivo model in rabbits. *Thrombosis and Haemostasis* 61(2), 307-313, 1989.

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XIV.LITERATUR

1. Elms, M.J. et al.: Rapid detection of cross-linked fibrin degradation products in plasma using monoclonal antibody-coated latex particles. J. Clin. Path. 85, 360-364, 1986.
2. Ballegeer, V. et al.: Fibrinolytic response to venous occlusion and fibrin fragment D-dimer levels in normal and complicated pregnancy. Thromb. Haemostas. 58, 1030-1032, 1987.
3. Declerck, P.V. et al.: Fibrinolytic response and fibrin fragment D-dimer in patients with deep vein thrombosis. Thromb. Haemostas. 58, 1025-1029, 1987.
4. Bounameaux, H. et al.: Measurement of plasma D-dimer for diagnosis of deep venous thrombosis. Am. J Clin. Path. 91, 82-85, 1989.
5. Hansson P.O. et al.: Can laboratory testing improve screening strategies for deep vein thrombosis at an emergency unit? J. Intern. Med.235, 143-151, 1994.
6. Holovet P et al.: Binding properties of monoclonal antibodies against human fragment D-dimer of cross-linked fibrin to human plasma clots in an in vivo model in rabbits. Thrombosis and Haemostasis 61(2), 307-313, 1989.

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XIV. BIBLIOGRAFIA

1. Elms, M.J. et al.: Rapid detection of cross-linked fibrin degradation products in plasma using monoclonal antibody-coated latex particles. *J. Clin. Path.* 85, 360-364, 1986.
2. Ballegeer, V. et al.: Fibrinolytic response to venous occlusion and fibrin fragment D-dimer levels in normal and complicated pregnancy. *Thromb. Haemostas.* 58, 1030-1032, 1987.
3. Declerck, P.V. et al.: Fibrinolytic response and fibrin fragment D-dimer in patients with deep vein thrombosis. *Thromb. Haemostas.* 58, 1025-1029, 1987.
4. Bounameaux, H. et al.: Measurement of plasma D-dimer for diagnosis of deep venous thrombosis. *Am. J Clin. Path.* 91, 82-85, 1989.
5. Hansson P.O. et al.: Can laboratory testing improve screening strategies for deep vein thrombosis at an emergency unit? *J. Intern. Med.* 235, 143-151, 1994.
6. Holovet P et al.: Binding properties of monoclonal antibodies against human fragment D-dimer of cross-linked fibrin to human plasma clots in an in vivo model in rabbits. *Thrombosis and Haemostasis* 61(2), 307-313, 1989.

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1. Elms, M.J. et al.: Rapid detection of cross-linked fibrin degradation products in plasma using monoclonal antibody-coated latex particles. *J. Clin. Path.* 85, 360-364, 1986.
2. Ballegeer, V. et al.: Fibrinolytic response to venous occlusion and fibrin fragment D-dímer levels in normal and complicated pregnancy. *Thromb. Haemostas.* 58, 1030-1032, 1987.
3. Declerck, P.V. et al.: Fibrinolytic response and fibrin fragment D dímer in patients with deep vein thrombosis. *Thromb. Haemostas.* 58, 1025-1029, 1987.
4. Bounameaux, H. et al.: Measurement of plasma D-dímer for diagnosis of deep venous thrombosis. *Am. J Clin. Path.* 91, 82-85, 1989.
5. Hansson P.O. et al.: Can laboratory testing improve screening strategies for deep vein thrombosis at an emergency unit? *J. Intern. Med.* 235, 143-151, 1994.
6. Holvet P et al.: Binding properties of monoclonal antibodies against human fragment D-dímer of cross-linked fibrin to human plasma clots in an in vivo model in rabbits. *Thrombosis and Haemostasis* 61(2), 307-313, 1989.

Helena Biosciences Europe
Queensway South
Team Valley Trading Estate
Gateshead
Tyne and Wear
NE11 0SD

Tel: +44 (0) 191 482 8440
Fax: +44 (0) 191 482 8442
Email: info@helena-biosciences.com
www.helena-biosciences.com

HL-2-0554P 2009/01 (9)