Protein S Antigen Rocket EID

Instructions For Use

REF 5359

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INTENDED PURPOSE

The Protein S procedure is used to determine the concentration of Protein S in human plasma as a free protein and as a complex with complement component C-4 binding protein (C4bp) using Laurell rocket electrophoresis^{1,2}.

Protein S is a vitamin K dependent plasma protein that acts as a potent cofactor for activated Protein C. In the coagulation process, it inhibits thrombin formation forming Factors V and VIII:C. Abnormally low Protein S levels are found in congenital Protein S deficiency^{3,4}, liver disorders and oral anticoagulant patients. The complexity of Protein S is illustrated by the formation of a dynamic equilibrium with C4b - binding protein. C4b - BP can be found in plasma in two different forms^{4,5}. The first is a low molecular weight form which exhibits no affinity for Protein S. The second is a high molecular weight form which can bind Protein S in a 1:1 complex.

Free Protein S represents about 40% (10 mg/l) of total Protein S. It acts as the carrier for anticoagulant activity and as the cofactor for activated Protein C. The C4b - BP bound Protein S normally represents about 60% (15 mg/l) of total Protein S and does not have any anticoagulant activity.

The Helena Protein S Rocket EID Procedure is performed in a 1% agarose gel medium containing an antiserum specific for Protein S. After the plasma specimens are applied to the wells in the agarose, electrophoresis is used to migrate the proteins into the antibody field. A rocket-shaped precipitin pattern forms along the axis of migration. The length of this rocket pattern is proportional to the antigen concentration.

Studies done by Dahlback⁷ showed that the addition of PEG to the plasma would precipitate the C4b - BP complex but not the free Protein S. The Helena Free Protein S Reagent resolves the free from the bound Protein S and thus allows for correlation between functional Protein S activity and immunologic levels of free Protein S.

WARNINGS AND PRECAUTIONS

The reagents contained in this kit are for *in-vitro* diagnostic use only - **DO NOT INGEST**. Wear gloves when handling all kit components. Refer to the product safety data sheets for risk and safety phrases and disposal information. Plasma products have been screened and found negative (unless otherwise stated on the kit box or vial) for the presence of Hepatitis B Antigen (HbsAg) HIV 1 and 2 antibody and HCV antibody, however they should be handled with the same precautions as a human plasma sample.

COMPOSITION

1. Protein S Antigen Rocket Plates (Cat. No. 5359)

Contains sheep or goat antibody to human Protein S incorporated into agarose in Tris-tricine buffer and sodium azide as preservative. To prevent the formation of toxic vapors, sodium azide should not be mixed with acidic solutions.

Preparation: Remove the plate from the protective packaging and allow 5-10 minutes for the agarose to reach 15...30°C.

2. Tris-tricine Buffer (Cat. No. 5358 - not included)

When diluted, the buffer contains 0.081M Tris and 0.024M tricine.

Preparation: Dilute one package of buffer to 1000ml with deionized water. The buffer is ready for use when all material is completely dissolved.

3. Rocket Stain (Cat. No. 5362 - not included)

Rocket Stain is Coomassie Brilliant Blue stain.

Preparation: Dissolve the contents of the vial in 450 ml deionized water, 450 ml methanol and 100 ml acetic acid. Mix thoroughly and filter before use if necessary.

4. Free Protein S Reagent (Cat. No. 5363)

The reagent contains 25% Polyethylene glycol. The reagent is ready for use as packaged.

5. Other Kit Components

Each reagent contains Instructions for Use, a rocket ruler and report form.

STORAGE AND SHELF-LIFE

Unopened reagents are stable until the given expiry date when stored under conditions indicated on the vial or kit label.

1. Protein S Antigen Rocket Plates (Cat. No. 5359)

Rocket Plates MUST be stored at 2...6°C and maintained in moist condition within the bag. **DO NOT FREEZE.** The plates are stable until the expiry date indicated on the package.

Signs of Deterioration: Discard the plate if dry in appearance or if the wells are not round. A crystalline appearance indicates the agarose has been frozen.

2. Tris-tricine Buffer (Cat. No. 5358)

The buffer should be stored at 15...30°C. Diluted buffer is stable for 2 months when stored at 15...30°C.

Signs of Deterioration: Discard packaged buffer if the material shows signs of dampness or discoloration. Discard diluted buffer if it becomes turbid.

3. Rocket Stain (Cat. No. 5362)

The stain should be stored at 15...30°C.

Signs of Deterioration: If methanol evaporation occurs, a metallic sheen will be visible on the stain surface. Discard the stain if it does not adequately stain protein rockets as described in this procedure.

4. Free Protein S Reagent (Cat. No. 5363)

The reagent should be stored at 2...6°C and is stable until the expiry date indicated on the label.

Signs of Deterioration: Discard the reagent if it shows any signs of microbial contamination.

ITEMS REQUIRED BUT NOT PROVIDED

Cat. No. 5358 Tris - Tricine Buffer

Cat. No. 5362 Rocket Stain

Cat. No. 5363 Free Protein S Reagent

Cat. No. 5185 SARP reference plasma

Cat. No. 4063 Titan Gel Chamber

Cat. No. 5014 Development Weight

Cat. No. 1520 EWS Power Supply

Cat. No. 5037 Blotter Pads

Cat. No. 9025 IEP VuBox

Cat. No. 5116 I.O.D. - Incubator, Oven, Drier

Cat. No. 1558 TITAN GEL Multi-staining Set

Cat. 6210, 6211 Microdispenser and Tubes (10ml)

Lint free tissue

Destain solution: Mix 200ml deionized water, 200ml methanol and 50ml glacial acetic acid.

SAMPLE COLLECTION AND PREPARATION

Plastic or siliconised glass should be used throughout. Blood (9 parts) should be collected into 3.2% or 3.8% sodium citrate anticoagulant (1 part). Separate plasma after centrifugation at 2000-3000xg for 15 minutes. Plasma should be kept at 2...6°C. Testing should be completed within 2 hours of sample collection, or plasma can be stored frozen at -20°C for 2 weeks or -70°C for one month. Thaw quickly at 37°C prior to testing. DO NOT keep at 37°C for more than 5 minutes.

STEP-BY-STEP PROCEDURE

- I. For step I follow the procedure according to Free Protein S Determination or Total Protein S Determination..

Free Protein S Determination

- i. All samples, including standard material and controls, must be pretreated by adding 15ml of Free Protein S Reagent to 85ml of plasma. These volumes may be altered, but the ratio should remain the same. DO NOT dilute the patient and control samples before pretreatment.
- ii. Vortex each specimen and place on ice for 30 minutes.
- iii. Centrifuge each specimen at a minimum of 1000 X G for 10 minutes. Reserve the supernatant for standard dilutions and patient dilutions.
- iv. To make the standard curve, dilute the standard material supernatant with 0.85% saline according to the following instructions.

Nominal Percent Activity	Dilution	Parts S.A.R.P.	Parts 0.85% Saline
100%	Use reconstituted S.A.R.P. undiluted		
50%	1:2	1	1
25%	1:4	1	3

- v. Go to Step 2

Total Protein S Determination

- i. Reconstitute one vial of S.A.R.P. with 1.0ml deionized water. Make dilutions for preparation of the Standard Curve as follows:

Nominal Percent Activity	Dilution	Parts S.A.R.P.	Parts 0.85% Saline
50%	1:2	1	1
25%	1:4	1	3
12.5%	1:8	1	7

- ii. Dilute each patient sample with 0.85% saline. Prepare a 1:2 dilution (1 part patient sample and 1 part saline) and a 1:4 dilution (1 part patient plasma and 3 parts saline). Additional dilutions may be necessary depending on the patient history. Suspected abnormal samples may need to be tested undiluted.
- iii. Go to Step 2

NOTE: Remove the plate from the refrigerator, allow approximately 5-20 minutes for the plate to reach 15...30°C and for excess moisture to be absorbed before use.

2. Pour 65ml of Tris-tricine buffer into each inner section of the chamber.
3. Remove any excess buffer from the wells of the plate. Excess moisture on the plate can result in poor rockets.

4. Apply 10ml of each patient sample or dilution to the appropriate well. Standard curve samples must be run on each plate. **NOTE:** Duplicate applications of patient samples are advisable. When applying the samples to the plate wells, do not allow the pipette tip to touch the sides of the wells as this may cause damage.
5. Allow 5 minutes for specimens to diffuse into the agarose.
6. Place the gel into the inner section of the chamber, agarose side down, by gently squeezing the gel into place. Position the gel(s) so that the edges of the agar are in the buffer and the wells are toward the cathode (-) side of the chamber.
7. Electrophorese the plates at a constant current of 16 mA per plate for 3 hours.
8. Following electrophoresis, remove the plates from the chamber. Discard the chamber buffer after each run. **NOTE:** Use the TITAN GEL Multi-Staining Set (Cat. No. 1558) as a staining and rinsing chamber.
9. Rinse the plate with deionized water and wash it in 0.85% saline overnight.
10. After the overnight wash, rinse the plate with deionized water.
11. Place the plate on a flat surface, agarose side up. Cover the agarose with a single, lint-free tissue.
12. Place 2-3 Blotter Pads and a Development Weight on the plate for 15 minutes and remove.
13. Dry the plate in a drying oven at 60...70°C for 10-20 minutes. DO NOT over dry plates. The plate will be transparent when completely dry. **NOTE:** If a dryer/oven is not available, the plates may be covered with wet lint-free tissues and allowed to dry at 15...30°C overnight or under a fan for 3 hours at 15...30°C as climate requires.
14. Following drying, stain the plate by immersing it in Rocket Stain for 20 minutes.
15. Place the plate in destain solution for 1-3 minutes. **NOTE:** Destaining is complete as soon as the background sufficiently clears in order to easily distinguish rocket peaks. Excessive destaining may fade the rockets making correct measurements difficult. If over destaining does occur, repeat Step 14 and stain the rockets again.
16. Rinse the plate twice in deionized water for 5-10 minutes each rinse.
17. Dry the plates at 37°C for 5 minutes or at 15...30°C until dry.
18. Place the plate on the Helena I.E.P.VuBox using a piece of white paper in the bottom of the VuBox for easier viewing of the rockets. Mark the apex of each rocket peak with a marker.
19. Using the Helena Rocket Ruler, measure the length of each peak in millimeters. The peak is measured from the top of each well to the apex of the rocket.
20. Plot the values of the standard curve versus each rocket height on the Rocket Antigen Report Form or on 3 cycle semi-logarithmic paper. Draw the line of best fit for the standard points. See Figures 1 and 2 in INTERPRETATION OF RESULTS for an example of a completed Rocket Plate and a standard curve drawn on a Rocket Antigen Report Form.

INTERPRETATION OF RESULTS

Read the patient values from the standard curve and multiply each by the appropriate dilution factor. If S.A.R.P. is used to prepare the standard curve, the patient value read from that curve must be multiplied by the assigned Protein Antigen value of the appropriate lot number of S.A.R.P. as well as the dilution factor.

Patient value from curve = 30%

Dilution factor = 2

S.A.R.P. assigned value = 98%

Actual Patient Factor

Protein S Antigen = $30\% \times 2 \times 0.98 = 58.8\%$

Patient samples with Protein S Antigen levels greater than the range of the standard curve, must be reassayed using the appropriate dilutions

Figure 1: Rocket patterns on a Protein S Antigen Rocket Plate. The lengths of the rockets (in millimeters) of the standard dilutions are used to prepare the standard curve. Patient results are read from the curve.



Figure 2: Representative standard curve prepared with S.A.R.P. on 3 cycle semi-logarithmic paper.

All decreases in Protein S, whether congenital or acquired, cause an increase in the risk of thromboembolism. Recurrent thrombotic episodes occur due to a decrease of anticoagulant potential³⁻⁶.

Protein S levels are decreased in:

- * Protein S congenital deficiency³⁻⁶.
 - type I, indicated by a decrease in the anticoagulant activity (free Protein S) and a normal to slightly decreased antigen level.
 - type II, indicated by a significant decrease in both the activity and antigen levels of Protein S.
- * Hepatic disorders.
- * Oral anticoagulant therapy.

QUALITY CONTROL

Each laboratory should establish a quality control program. Normal and abnormal control plasmas should be tested prior to each batch of patient samples, to ensure satisfactory instrument and operator performance. If controls do not perform as expected, patient results should be considered invalid.

Helena BioSciences supply the following controls available for use with this product:

Cat. No. 5301 SAC 1

Cat. No. 5302 SAC 2

REFERENCE VALUES

Reference values can vary between laboratories depending on the techniques and systems in use. For this reason each laboratory should establish its own normal range.

Protein S antigen values are usually expressed in relative percentages compared to a pooled normal plasma standard.

Total Proteins: The normal range study consisted of fifty (50) normal healthy adults, 25 males & 25 females. There has been no indication of any difference in Protein S antigen range levels between healthy males & females. The expected normal range for Protein S antigen run by rocket electrophoresis was 60-150%.

These specimens were from donors who are aged between 18 and 65 years old, in good health; no surgery within 6 months, have normal haemoglobin, haematocrit and temperature and are taking no medication, i.e., antibiotics.

Free Protein S: Normal free Protein S levels established by this method have been reported from 57% - 120%⁵. Helena tested 31 plasmas and determined an actual range of specimen concentration from 42% to 147%. Each laboratory should establish its own normal for this procedure but be aware that there is currently no accepted international standard for Protein S antigen.

LIMITATIONS

Heparin samples with levels up to 10 IU/ml do not affect the results of the Helena Protein S antigen assay. Patients on oral anticoagulant therapy do show marked decreases in the protein S antigen level. The free Protein S level does not affect the measurement of total Protein S. No other substances have been found to interfere with this procedure.

PERFORMANCE CHARACTERISTICS

Helena BioSciences or their representatives have determined the following performance characteristics as a guideline. Each laboratory should establish its own performance data.

Precision

Studies were done with a control material used in replicate on one plate.

	Within Run:		Between Run:	
	Total S	Free S	Total S	Free S
\bar{x}	98.1	50.8	98.4	52.2
SD	2.8	2.7	2.8	3.5
CV%	2.3	5.2	5.6	6.7
n	8	5	18	15

Comparison

Correlation studies were done using the Helena Protein S system and a reference method. The specimens used met the same criteria as those for the normal range study. The following data was derived.

	Total S		Free S
n	= 49	31	Y = Helena Protein S
Y	= 0.917X + 6.975	0.981X + 0.498	X = Reference Method
r	= 0.849		0.906

Sensitivity

Sensitivity for this procedure is expressed in relative percentages compared to a pool of normal plasma, such as S.A.R.P. Patient values greater than the highest value on the normal standard curve must be diluted and retested. The sensitivity threshold for the Helena Protein S Kit is 6% of the normal human level.

Specificity

Specificity was done on known Protein S deficient plasmas and verified using a Double Diffusion (Ouchterlony) plate and various dilutions of antisera associated with coagulation proteins. The Helena Protein S Kit utilizes goat antihuman Protein S in optimum concentrations.

BIBLIOGRAPHY

1. Laurell, C.B., Electroimmuno Assay, Scan J Clin Lab Invest, 29 Suppl 124:21-37, 1972.
2. Laurell, C.B., Quantitative Estimation of Proteins by Electrophoresis in Agarose Gel Containing Antibodies. Anal Biochem, 15:45-52, 1966.
3. Schwarz H.P., et al, Plasma protein S deficiency in familial thrombotic disease. Blood, 64:1297-1300, 1984.
4. Comp P.C., et al., Familial protein S deficiency is associated with recurrent thrombosis. J Clin Invest, 74:2082-2088, 1984.
5. Comp P.C., et al., An abnormal plasma distribution of protein S occurs in functional protein S deficiency. Blood, 67:504-508, 1986.
6. Kamiya T., et al., Inherited deficiency of protein S in a Japanese family with recurrent venous thrombosis: a study of three generations. Blood, 67:406-410, 1986.
7. Dahlback, B., Purification of human vitamin K-dependent protein S and its limited proteolysis by thrombin. Biochem Jour, 209:827, 1983.

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BIBLIOGRAPHIE

1. Laurell, C. B., Electroimmuno Assay, Scan J Clin Lab Invest 29 Suppl 124 : 21-37, 1972.
2. Laurell, C. B., Quantitative Estimation of Proteins by Electrophoresis in Agarose Gel Containing Antibodies. Anal Biochem, 15 : 45-52, 1966.
3. Schwarz H. P., et al, Plasma protein S deficiency in familial thrombotic disease. Blood, 64 : 1297-1300, 1984.
4. Comp P. C., et al., Familial protein S deficiency is associated with recurrent thrombosis. J Clin Invest, 74 : 2082-2088, 1984.
5. Comp P. C., et al., An abnormal plasma distribution of protein S occurs in functional protein S deficiency. Blood, 67 : 504-508, 1986.
6. Kamiya T., et al., Inherited deficiency of protein S in a Japanese family with recurrent venous thrombosis: a study of three generations. Blood, 67 : 406-410, 1986.
7. Dahlback, B., Purification of human vitamin K-dependent protein S and its limited proteolysis by thrombin. Biochem Jour, 209 : 827, 1983.

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