en

INTENDED PURPOSE

For use as a control in the Ristocetin CoFactor Assay of human plasma.

WARNINGS AND PRECAUTIONS

The reagents contained in this kit are for in-vitro diagnostic use only - **DO NOT INGEST**. Wear gloves when handling all kit components. Refer to the product safety data sheets for risk and safety phrases and disposal information. Plasma products have been screened and found negative (unless otherwise stated on the kit box or vial) for the presence of Hepatitis B Antigen (HbsAg) HIV 1 and 2 antibody and HCV antibody; however they should be handled with the same precautions as a human plasma sample.

COMPOSITION

1. Ristocetin CoFactor Abnormal Control Plasma (5 × 0.5ml)

Each vial contains a lyophilised and stabilised preparation of citrated human plasma with a Ristocetin CoFactor level of less than 50%.

Preparation: Reconstitute each vial with 0.5ml of distilled or deionised water, swirl gently and allow to stand for 10 minutes at room temperature. Mix well before use (do not shake).

2. Other kit components

Each kit contains Instructions For Use and lot specific reference values insert.

STORAGE AND SHELF-LIFE

1. Ristocetin CoFactor Abnormal Control Plasma

Unopened vials are stable until the given expiry date when stored under conditions indicated on the vial label. Avoid contamination. Once reconstituted, the plasma is stable for 8 hours at 2...6°C or 1 month at -20°C.

STEP BY STEP PROCEDURE

Each control should be treated in the same manner as the unknown specimen in accordance with the instructions outlined in each particular test protocol.

REFERENCE VALUES

Reference values can vary between laboratories depending on the techniques and systems in use. For this reason each laboratory should establish its own reference ranges.

Lot specific expected values are provided with each pack of controls.

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