### INTENDED PURPOSE

Thrombin Reagent is intended for use in the quantitative determination of fibrinogen in human plasma using the Clauss Method.

Clauss' developed a simple method for the quantitative determination of fibrinogen by measuring the clotting time of dilute plasma after the addition of thrombin (>30 NIH units/ml). This clot time is proportional to the fibrinogen concentration. Levels of fibrinogen can increase as a result of inflammation, pregnancy or oral contraceptive use<sup>2</sup>. Decreased levels can be found in certain states such as liver disease and DIC. Congenital deficiencies include afibrinogenaemia (no detectable fibrinogen), hypofibrinogenaemia (<1 mg/ml) and dysfibrinogenaemia (abnormal fibrinogen molecule).

### WARNINGS AND PRECAUTIONS

The reagents contained in this kit are for in-vitro diagnostic use only - **DO NOT INGEST**. Wear gloves when handling all kit components. Refer to the product safety data sheets for risk and safety phrases and disposal information.

Plasma products have been screened and found negative (unless otherwise stated on the kit box or vial) for the presence of Hepatitis B Antigen (HbsAg) HIV 1 and 2 antibody and HCV antibody; however they should be handled with the same precautions as a human plasma sample.

### COMPOSITION

- Thrombin Reagent**  
**REF 5374 (10 x 2ml)**  
**REF 5378 (10 x 5ml)**  
Each vial contains approximately 100 NIH units/ml bovine thrombin with stabilisers.  
**Preparation:** Reconstitute each vial with 2.0ml (for Cat. No. 5374) or 5.0ml (for Cat. No. 5378) of purified water. Swirl gently and allow to stand for 15 minutes. Mix well before use. Do not shake.
- Other kit components**  
Each kit contains Instructions For Use.

### STORAGE AND SHELF-LIFE

Unopened vials are stable until the given expiry date when stored under conditions indicated on the vial or kit label.

- Thrombin Reagent**  
Once reconstituted, the reagent is stable for 8 hours at room temperature or 1 week at 2...6°C.

### ITEMS REQUIRED BUT NOT PROVIDED

Cat. No. 5379 Fibrinogen Calibrator 10 x 1ml or  
Cat. No. 5185 SARP 10 x 1ml  
Cat. No. 5375 Owren's Veronal Buffer 10 x 25ml

Helena BioSciences Europe also provides the complete Fibringoen Assay kit, Cat. No. 5376.

### SAMPLE COLLECTION AND PREPARATION

Plastic or siliconised glass should be used throughout. Blood (9 parts) should be collected into 3.2% or 3.8% sodium citrate anticoagulant (1 part). Separate plasma after centrifugation at 2000 - 3000 x g for 15 minutes. Plasma should be kept at 2...6°C. Testing should be completed within 2 hours of sample collection, or plasma can be stored frozen at -20°C for 2 weeks or -70°C for one month. Thaw quickly at 37°C prior to testing. Do not keep at 37°C for more than 5 minutes.

### STEP BY STEP PROCEDURE

Helena BioSciences Thrombin reagent is suitable for use on any coagulation analyser capable of performing Clauss Fibrinogen assays. Refer to the Instrument Operators Manual and/or application note for appropriate instructions.

Prepare all reagents as instructed under 'composition'..

- Standard Curve Preparation**  
Note: A new standard curve should be used when reagent lot numbers change or if expected values in quality control testing change.  
Prepare the following dilutions in Owren's Veronal Buffer:

Tube	Dilution	Fibrinogen Calibrator (ml)	Buffer (ml)
1	1+4	0.2	0.8
2	1+9	0.1	0.9
3	1+19	0.1	1.9
4	1+29	0.1	2.9
5	1+39	0.1	3.9

Mix without shaking.
- Patient Sample Preparation**  
Prepare 1+9 dilutions of the patient plasma or control plasma in Owren's Veronal Buffer. Mix without shaking.
- Testing**  
Perform all tests in duplicate.
  - Pipette 0.2ml of standard, patient or control dilution into a reaction tube and incubate at 37°C for 2 minutes.
  - Add 0.1ml of Thrombin Reagent 15...30°C.
  - Determine the clot time to the nearest 0.1 seconds.
  - Plot mean Standard Clot Time (Y-axis) versus Fibrinogen Level (X-axis) on the Fibrinogen graph paper.  
A straight line should be obtained. Assign the Fibrinogen Calibrator reference value to the 1+9 dilution to allow direct interpolation of patient and control values from the standard curve.

### QUALITY CONTROL

Each laboratory should establish a quality control program. Normal and abnormal control plasmas should be tested prior to each batch of patient samples, to ensure satisfactory instrument and operator performance. If controls do not perform as expected, patient results should be considered invalid.

Helena BioSciences supply the following controls available for use with this product:
REF 5301 Sac-1 Normal Control Plasma 10x1ml
REF 5302 Sac-2 Abnormal Control Plasma 10x1ml
REF 5186 Norm-Trol 1 Control Plasma 10x1ml
REF 5187 Ab-Trol 2 Control Plasma 10x1ml
REF 5183 Ab-Trol 3 Control Plasma 10x1ml

### REFERENCE VALUES

Reference values can vary between laboratories depending on the techniques and systems in use. For this reason each laboratory should establish it's own normal range. Expected values for fibrinogen in healthy adults are 150-350 mg/dl (1.5-3.5 g/L)<sup>3,4</sup>.

### LIMITATIONS

Heparin levels >0.6 units/ml and fibrinolytic degradation products >100mg/ml may cause falsely low fibrinogen quantitation. If values fall outside the standard curve values for the patient samples, re-assay using an appropriate dilution to bring values into the standard range.

### PERFORMANCE CHARACTERISTICS

Helena BioSciences or their representatives have determined the following performance characteristics as a guideline. Each laboratory should establish it's own performance data. The Fibrinogen assay is designed to give a linear calibration from 1.5 - 6.5g/L.

#### REPRODUCIBILITY:

Within assay reproducibility: Fibrinogen (g/L)=3.64 CV=4.87 n=10	
Between assay reproducibility: Fibrinogen (g/L)=3.61 CV=3.78 n=10	

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