INTENDED PURPOSE

The Helena BioSciences Thrombin Clotting Time Reagent contains bovine thrombin. The reagent can be used manually or on semi-automated and automated instruments. The test is commonly applied to detect various sources of interference with normal blood coagulation. Prolongation of the thrombin clotting time can be taken as a qualitative indication of abnormal fibrinogen levels (high or low), or the presence of interfering substances such as FDP's or heparin. Quantitative evaluation of the possible causes of prolonged thrombin clotting time should be performed as follow-up studies, such as aPTT or chromogenic assay for heparin, Clauss fibrinogen, FDP determinations, heparin neutralisation by protamine sulphate or polybrene, normal plasma mixing studies' or reptilase assay' to distinguish between hypofibrinogenaemia and FDP effects.

WARNINGS AND PRECAUTIONS

The reagents contained in this kit are for in-vitro diagnostic use only - **DO NOT INGEST**. Wear gloves when handling all kit components. Refer to the product safety data sheets for risk and safety phrases and disposal information. Plasma products have been screened and found negative (unless otherwise stated on the kit box or vial) for the presence of Hepatitis B Antigen (HbsAg) HIV 1 and 2 antibody and HCV antibody; however they should be handled with the same precautions as a human plasma sample.

COMPOSITION

- Thrombin Clotting Time Reagent**
REF 5392 (10 x 2ml)
REF 5377 (10 x 5ml)
 Each vial contains a lyophilised preparation of bovine thrombin with buffers and stabilisers. The reconstituted reagent contains approximately 3NIH units/ml of thrombin. The reagent should be a white lyophilised plug.
Preparation: Reconstitute the reagent with the recommended volume of purified water (5ml for REF 5377 or 2ml for REF 5392). Allow to stand for 5 minutes then mix gently by inversion and transfer to a plastic tube. Reconstituted reagent should be a clear, colourless solution.
- Other kit components**
 Each kit contains Instructions For Use.

STORAGE AND SHELF LIFE

Unopened reagents are stable until the given expiry date when stored under conditions indicated on the vial or kit label.

- Thrombin Clotting Time Reagent**
 Reconstituted reagent should be stored at 2...6°C and is stable for 8 hours, AC-4 Pos 27 @ 37°C = 6hrs.

SAMPLE COLLECTION AND PREPARATION

Plastic or siliconised glass should be used throughout. Blood (9 parts) should be collected into 3.2% or 3.8% sodium citrate anticoagulant (1 part). Separate plasma after centrifugation at 2000 - 3000 x g for 15 minutes. Plasma should be kept at 2...6°C. Testing should be completed within 2 hours of sample collection, or plasma can be stored frozen at -20°C for 2 weeks or -70°C for one month. Thaw quickly at 37°C prior to testing. Do not keep at 37°C for more than 5 minutes.

STEP-BY-STEP PROCEDURE

Any high quality electro-mechanical or photo-optical coagulation instrument capable of performing Thrombin Clotting Time testing may be used.

- Collect and prepare the blood sample according to the directions outlined in **SAMPLE COLLECTION AND PREPARATION**.
- Reconstitute the control plasma according to the package insert provided with each control.
- Prepare the reagent for use in the procedure according to the reconstitution instructions in the **COMPOSITION** section.
- Perform all tests in duplicate. Calculate the mean clotting time of the duplicate determinations to the nearest 0.1 seconds. Individual values should be within ±5% of the mean value.

- Automated Methods, AC-4. (REF 5392 x400 samples, REF 5377 x1000 samples)**
 Refer to the appropriate AC-4 Operators Manual for detailed instructions. Prepare all reagents as instructed under 'composition'.

START =5s	METHOD =Coag	MIX=0	PATIENT (Pat) vol=100µL Pos=35 (CP)	Def. (DP) Vol=0µL	Plasma
INCUB =0s	MATH =lin	CLEAN=0	BUFFER (BUF) Vol=0µL	R1 Vol=0µL	
RUNTIME =120s	SENS =2	T-CORR=0 S-CORR=0	CLEAR (CLR) Vol=0µL	R2 Vol=50µL Pos=27	

B. Manual Method

- Pipette 0,2ml of the patient plasma or control plasma into a reaction tube.
- Incubate at 37°C for 3 minutes.
- Pipette 0,1ml of Thrombin Clotting Time Reagent into the reaction tube containing patient or control plasma whilst simultaneously starting a timer.
- Record the time taken for clot formation to the nearest 0.1 second.

QUALITY CONTROL

Each laboratory should establish a quality control program. Normal and abnormal control plasmas should be tested prior to each batch of patient samples, to ensure satisfactory instrument and operator performance. If controls do not perform as expected, patient results should be considered invalid.

Helena BioSciences supply the following controls available for use with this product:

REF 5186	Norm-Trol 1	(10 x 1ml)
REF 5187	Ab-Trol 2	(10 x 1ml)
REF 5183	Ab-Trol 3	(10 x 1ml)

REFERENCE VALUES

Reference values can vary between laboratories depending on the techniques and systems in use. For this reason each laboratory should establish it's own normal range.

INTERPRETATION OF RESULTS

The results of the Thrombin Clotting Time test should be reported to the nearest 0.1 seconds. The normal range (usually the mean ±2 standard deviations) should be established by each laboratory. Results outside the normal range should be considered abnormal and follow-up testing performed.

LIMITATIONS

Expected values for the Thrombin Clotting Time test will vary from one laboratory to another depending on the technique used. The method of clot detection, temperature, pH, collection technique, type of anticoagulant and time and method of sample storage are all very important. Plasma sample collection and storage conditions should be standardised and carefully controlled. Unexpected results should be confirmed by additional testing. As well as the cause of elongated Thrombin Clotting Times indicated in the summary section, a recent report has suggested that many systemic amyloidosis patients with bleeding complications may have a circulating inhibitor which prolongs the Thrombin Clotting Time⁴. Also, therapeutic levels of heparin may entirely abolish clotting in the Thrombin Clotting Time test, although neutralisation with protamine sulphate or polybrene should correct the Thrombin Clotting Time¹.

PERFORMANCE CHARACTERISTICS

Helena BioSciences or their representatives have determined the following performance characteristics:

Within assay reproducibility		
Clot Time (s) CV(%)	n	
12.5	1.43	10
Between assay reproducibility		
Clot Time (s) CV(%)	n	
12.5	1.29	10

BIBLIOGRAPHY

- Laposata et al. The Clinical Haemostasis Handbook, Yearbook Medical Publishers Inc., p219, 1989.
- Thompson, A.R. and Harker, L.A. Manual of Haemostasis and Thrombosis. 3rd Ed., F.A. Davis Co., p62, 1983.
- DeMott, W.R., in: Laboratory Test Handbook, 2nd Ed., Jacobs D.S. et al Eds., Lexi-Comp Inc., p432-433, 1990.
- Gastineau, D.A. et al. 'Inhibitor of the Thrombin Time in Systemic Amyloidosis: A Common Coagulation Abnormality' Blood, 1991, 77: 2637-2640.

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