## INTENDED PURPOSE

Lupus Anticoagulants (LA's) are antibodies of the IgG and IgM type which are directed against a variety of anionic phospholipids. The presence of LA's in plasma is increasingly associated with a variety of haemostatic problems such as recurrent foetal loss, thrombocytopaenia, unexplained thrombosis and neurological disorders. LA's prolong phospholipid dependant in vitro clotting assays such as the APTT. The Helena BioSciences DRVVT-Screen kit is intended for the qualitative determination of Lupus Anticoagulants (LA's) in human plasma. Russell's Viper Venom directly activates Factor X to Factor Xa in the presence of phospholipid and calcium, leading to detectable clot formation in plasma. The DRVVT-Screen kit is more sensitive for LA's than the APTT. The DRVVT-Screen kit is intended to be used in conjunction with the DRVVT-Confirm kit (REF 5485).

## WARNINGS AND PRECAUTIONS

The reagents contained in this kit are for in-vitro diagnostic use only - **DO NOT INGEST**. Wear gloves when handling all kit components. Refer to the product safety data sheets for risk and safety phrases and disposal information.

## COMPOSITION

### 1. DRVVT Screen Reagent (10 x 2ml)

Each vial contains a proprietary mixture of Russell's Viper Venom co-lyophilised with calcium chloride and phospholipid. Preparation: Reconstitute each vial with 2ml of distilled/deionised water. Allow to stand for 10 minutes and mix well before use (Do not shake).

### 2. Other kit components

Each kit contains Instructions For Use .

## STORAGE AND SHELF-LIFE

### 1. DRVVT Screen Reagent

Unopened vials are stable until the given expiry date when stored under conditions indicated on the vial or kit label. Reconstituted vials are stable for 24 hours at 15...30°C, 5 days at 2...6°C or 2 weeks at -20°C, AC-4 Pos 37 @ 37°C = 24hrs. . The reagent should be frozen in plastic test tubes and thawed at 37°C before use.

## ITEMS REQUIRED BUT NOT PROVIDED

REF 5485 DRVVT-Confirm Kit 10 x 1ml  
REF 5499 Normal Plasma (e.g. Norm-Trol 1) 10 x 3ml  
REF 5186 Normal Plasma (e.g. Norm-Trol 1) 10 x 1ml

## SAMPLE COLLECTION AND PREPARATION

Plastic or siliconised glass should be used throughout. Blood (9 parts) should be collected into 3.2% or 3.8% sodium citrate anticoagulant (1 part). Separate plasma after centrifugation at 2000xg - 3000xg for 15 minutes. Residual platelets should be removed by either filtration through a 0.22mm disposable filter, or re-centrifugation. Plasma should be kept at 2...6°C. Testing should be completed within 2 hours of sample collection, or plasma can be stored frozen at -20°C for 2 weeks or -70°C for one month. Thaw quickly at 37°C prior to testing. Do not keep at 37°C for more than 5 minutes.

## STEP BY STEP PROCEDURE

### 1. Automated Methods, AC-4.(x260 samples)

Refer to the appropriate AC-4 Operators Manual for detailed instructions.

START =11s	METHOD =Coag	MIX=0	PATIENT (Pat) vol=75µL pos=36 (CP)	Def. (DP) Vol=0µL	Plasma
INCUB =0s	MATH =lin	CLEAN=0	BUFFER (BUF) Vol=0µL	R1 Vol=0µL	
RUNTIME =240s	SENS =1	T-CORR=0 S-CORR=0	CLEAR (CLR) Vol=0µL	R2 Vol=75µL Pos=37(DRVVT -Screen)	

### 2. Manual Method

- Pre-warm sufficient reconstituted reagent to 37°C.
- Pipette 0.2ml of patient or control plasma into a reaction tube. Incubate at 37°C for 2 minutes.
- Add 0.2ml of pre-warmed DRVVT-Screen Reagent and start a timer.
- Measure the clot formation time to the nearest 0.1 seconds.
- Calculate the ratio

Patient DRVVT-Screen Clot Time  
Norm-Trol 1 DRVVT-Screen Clot Time

### 3. Automated Methods

The DRVVT-Screen test may be performed on most automated instruments. The instrument should be programmed to deliver equal volumes of plasma and reagent and both sample and reagent should be pre-warmed to 37°C for at least 2 minutes prior to mixing. Reagent tubing and reservoirs should be thoroughly rinsed and cleaned prior to performing other clotting or chromogenic tests. Refer to the appropriate Operators Manual for specific detailed instructions.

## INTERPRETATION OF RESULTS

Results are best expressed as a ratio relative to the clot times obtained on Norm-Trol 1. Both DRVVT-Screen and DRVVT-Confirm results can be 'normalised' in this way, reducing the effects of instrument variability and reagent batch differences. Results of mixing tests can be treated in the same way.

If the clot time of the patient samples are greater than 3 standard deviations above the mean of the normal range, a lupus anticoagulant may be present. In this case, the plasma should be re-tested after mixing 1:1 with Norm-Trol 1 as well as testing with the DRVVT-Confirm kit (Cat. No. 5485). If the DRVVT-Screen clotting time of the patient plasma mixed 1:1 with Norm-Trol 1 is still greater than 3 standard deviations above the mean of the normal range, a lupus anticoagulant may be present.

If the DRVVT-Screen clotting time of the patient plasma mixed 1:1 with Norm-Trol 1 is corrected to within the normal range, a factor deficiency (II, V or X) is most likely.

The Scientific and Standardisation Sub-Committee for the Standardisation of Lupus Anticoagulants of the International Society of Thrombosis and Haemostasis has recommended that the diagnosis of lupus anticoagulant be made when the DRVVT of a test plasma mixed with normal plasma is greater than 3 standard deviations from the mean normal (non-LA) plasma DRVVT time.

The use of the DRVVT-Confirm kit allows discrimination between LA, factor deficiency and other inhibitors.

## REFERENCE VALUES

Reference values can vary between laboratories depending on the techniques and systems in use. For this reason each laboratory should establish it's own normal range (range = mean ± 3sd's) using a minimum of twenty (20) healthy, male and female donors spanning the adult age range. If only frozen patient samples are tested, the normal reference range should be determined on frozen normal samples.

The normal reference range (mean ± 3sd's) determined at Helena BioSciences for the DRVVT-Screen test was 33.7 ± 8.1 seconds (range 25.6 - 41.8 seconds).

## LIMITATIONS

Plasma deficiencies of Factors II, V or X may lead to abnormal results in neat plasma. Mixing studies should correct this. Plasma from patients with the following may give abnormal results when the plasma is tested neat, and these samples may not correct in mixing studies: heparin (>1U/ml), oral anticoagulants, disseminated intravascular coagulation (DIC). Care must be taken to remove residual platelets from plasma by filtration or centrifugation, as platelet derived phospholipid can interfere with the test.

## QUALITY CONTROL

Each laboratory should establish a quality control program. Normal and abnormal control plasmas should be tested prior to each batch of patient samples, to ensure satisfactory instrument and operator performance. If controls do not perform as expected, patient results should be considered invalid.

Helena BioSciences supply the following controls available for use with this product:

REF 5486 DRVVT Positive Control 1 x 1ml

## PERFORMANCE CHARACTERISTICS

Helena bioSciences or their representatives have determined the following performance characteristics as a guideline. Each laboratory should establish it's own performance data. Within run and between run CV's are expected to be <5%.

## BIBLIOGRAPHY

- National Committee for the National Laboratory (NCCLS) Standards: Collection transport and preparation of blood specimens for coagulation testing and performance of coagulation assays. Document H21-A2, vol. 11, No. 23, 1991.

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