## INTENDED PURPOSE

The first standardized one-stage prothrombin time test was developed by Dr. Armand Quick in 1935<sup>2</sup>. It is used for the induction and monitoring of oral anticoagulant therapy<sup>2,3</sup> and can be used to assess the protein synthesis capability of the liver in chronic or acute hepatic disorders. The Helena BioSciences Capillary Reagent is of rabbit brain origin but resembles human BCT in its low International Sensitivity Index (ISI). The ISI of Capillary Reagent is approximately 1.1 (using manual calibration techniques) and is calibrated against the WHO international reference preparation<sup>1</sup>. Capillary Reagent is particularly suited to the monitoring of oral anticoagulant therapy.

Tissue thromboplastin, in the presence of calcium ions, is an activator which initiates the extrinsic pathway of coagulation. When a mixture of tissue thromboplastin and calcium ions is added to capillary or citrated blood, the clotting mechanism is activated, leading to a fibrin clot. If a deficiency exists within the extrinsic pathway, the time required for clot formation will be prolonged depending on the severity of the deficiency.

## WARNINGS AND PRECAUTIONS

The reagents contained in this kit are for In-Vitro diagnostic use only - **DO NOT INGEST**. Wear gloves when handling all kit components. Refer to the product safety data sheets for risk and safety phrases and disposal information.

## COMPOSITION

- Capillary Reagent (10 x 2.5ml)** Each vial contains rabbit brain thromboplastin suspension. Preparation: The contents of the vial should be mixed well before and during use.
- Calcium Chloride (10 x 2.5ml)** Each vial contains 0.010M calcium chloride solution.
- Other Kit Components** Each kit Contains Instructions For Use and lot specific reference values insert.

## STORAGE AND SHELF LIFE

Unopened reagents are stable until the given expiry date when stored under conditions indicated on the vial or kit label.

- Capillary Reagent** **Opened vials:** 1 month at 2...6°C and 5 days at 15...30°C. The working reagent (mixed with calcium chloride) is stable for 10 days at 2...6°C or 5 days at 15...30°C. **DO NOT FREEZE.** Signs of **Deterioration:** The reagent is a fine suspension of rabbit brain tissue. Large clumps of particles or changes in expected values may indicate product deterioration.
- Calcium Chloride** **Opened vials:** 1 month at 2...6°C or 5 days at 15...30°C. Signs of **Deterioration:** The reagent is a clear, colourless solution. Particulate contamination, cloudiness or a change in expected performance may indicate product deterioration.

## ITEMS REQUIRED BUT NOT PROVIDED

REF 90.200.0001 CoaData 501WB Fibrin Timer or any other whole blood coagulometer

## SAMPLE COLLECTION AND PREPARATION

Throughout the procedure, all test tubes, syringes and pipettes must be plastic or siliconized glass. Venous blood should be collected into 3.2% sodium citrate – 9 parts blood + 1 part sodium citrate. Capillary blood should be collected into the first drop of blood produced from a skin puncture of suitable depth. Testing on venous blood should be completed within 4 hours of sample collection and the sample kept at room temperature<sup>2</sup>.

## STEP-BY-STEP PROCEDURE

For accurate INR reporting, it is recommended to determine the laboratory specific ISI of the reagent with the testing system in use. The Helena BioSciences Calibrant Plasma Set (REF 5519) is recommended for this purpose<sup>3</sup>. This should be performed for each new reagent batch. The Helena BioSciences INR Reference Set (REF 5490) should be used to check for shifts in the local system ISI which have been noted with changes in laboratory temperature and post instrument servicing, amongst other local variances.

### Running instructions for the CoaData 501 WB:

- Prior to use the instrument should be allowed to equilibrate to 37°C.
- Pre-warm sample cups (including str bars) in the appropriate positions.
- Mix reagent and calcium chloride 1 + 1 and pre-warm in the reagent well on the instrument.
- Once warm the instrument sit at running temperature, press ‘Esc’, the screen will show <1PT-WB>.
- Press ‘Enter’, the screen will show ‘Start’.
- Insert a pre-warmed cup (containing 150µl reagent) into the channel and press ‘Start’.
- Wait for incubation, set at 1 second.
- At the prompt, ‘Go-S’, add 25µl blood (15µl plasma) into the cup in the measuring channel and simultaneously press ‘Start’.
- Wait for clot time/INR.
- After recording results remove cup and press ‘Reset’.

### Running instructions for the Thrombostat:

- Prior to use the instrument should be allowed to equilibrate to 37°C.
- Mix Reagent and Calcium Chloride (1 + 1) and pre-warm in the incubation position for reagents for 5 minutes.
- i) Thrombostat 1 users:** Once temperature is reached the Thrombostat is ready to run a PT. **ii) Thrombostat 2 users:** Press ‘Test’ until ‘PT test is selected. Place cuvette(s) into the incubation position for sample(s) and dispense 1 ball bearing into each.
- Whole Blood Method:**
  - Pipette 150µl of pre-warmed Reagent and Calcium Chloride (1 + 1) mix into cuvette(s). Incubate for 2 minutes:
    - Thrombostat 1 users:** Press the ‘Reset/Incubation’ to start the timer.
    - Thrombostat 2 users:** Press cuvettes to activate timer, the Thrombostat will beep to indicate the end of incubation.
  - Following incubation transfer cuvette(s) to the measuring channel position.
- Thrombostat 1 users:** Press ‘Start’
  - Thrombostat 2 users:** Press ‘Reset/Start’ The instrument will count down to the start of measurement (Thrombostat 2 users: channel 1). Just before measurement starts pipette 25µl of whole blood into the cuvette. **Thrombostat 2 users:** To run a duplicate press ‘Reset/Start’ again to begin measurement in channel 2.

### Plasma Method:

- Pipette 50µl of plasma into cuvette(s). Incubate for 2 minutes:
  - Thrombostat 1 users:** Press the ‘Reset/Incubation’ to start the timer.
  - Thrombostat 2 users:** Press cuvettes to activate timer, the Thrombostat will beep to indicate the end of incubation.
- Following incubation transfer cuvette(s) to the measuring channel position.
- i) Thrombostat 1 users:** Press ‘Start’
  - Thrombostat 2 users:** Press ‘Reset/Start’ The instrument will count down to the start of measurement (Thrombostat 2 users: channel 1). Just before measurement starts pipette 250µl of Reagent and Calcium Chloride (1 + 1) mix into the cuvette. **Thrombostat 2 users:** To run a duplicate press ‘Reset/Start’ again to begin measurement in channel 2. **NOTE:** If using the start pipette the measuring time is automatically started in step (c) eliminating the need to press ‘Start’ or ‘Reset/Start’.
- Wait for clot time/INR
- i) Thrombostat 1 users:** Record results, remove cuvette and press ‘Reset/Incubation’ to erase the results of the former test. The instrument is now ready for next test. **ii) Thrombostat 2 users:** After results have been printed remove cuvettes and press ‘Reset/Start’ to erase the results of the former test. The instrument is now ready for next test.

## INTERPRETATION OF RESULTS

Results should be reported to the nearest 0.1 seconds and duplicates should agree within 5% of each other.

INR values can be calculated using the following formula:

INR = (PT Time Patient) / (Mean Normal PT Time)<sup>ISI</sup>

%PT values can be interpolated from the calibration graph (%PT of PT Calibration plasmas versus measured clot time) which should be joined point-to-point when plotted on log-log graph paper.

## QUALITY CONTROL

Each laboratory should establish a quality control program. Normal and abnormal control plasmas should be tested prior to each batch of patient samples, to ensure satisfactory instrument and operator performance. If controls do not perform as expected, patient results should be considered invalid.

Helena BioSciences supply the following controls available for use with this product:

REF 5186 / 5499 Norm-Trol 1  
REF 5187 Ab-Trol 2  
REF 5183 Ab-Trol 3

## REFERENCE VALUES.

Reference values can vary between laboratories depending on the techniques and systems in use. For this reason each laboratory should establish it's own normal range.

This is particularly important for local ISI calibration.

Using the CoaData 501WB, normal values ranging from 19.5 - 29.0 seconds are typical.

### Therapeutic Ranges

The British Society for Haematology<sup>9</sup> recommend the following INR ranges for various clinical conditions:

| INR     | Clinical State  |
|---------|---|
| 2.0-2.5 | Prophylaxis of deep vein thrombosis, including surgery on high risk patients (2.0-3.0 for hip surgery and fractured femur operations).  |
| 2.0-3.0 | Prevention of pulmonary embolism. Systemic embolism. Prevention of venous thrombo-embolism in myocardial infarction. Mitral stenosis with embolism. Transient ischaemic attacks. Atrial fibrillation. |
| 3.0-4.5 | Recurrent deep vein thrombosis and pulmonary embolism. Arterial disease, including myocardial infarction. Mechanical prosthetic heart valves.   |

## PERFORMANCE CHARACTERISTICS

Helena BioSciences or their representatives have determined the following performance characteristics as a guideline. Each laboratory should establish it's own performance data.

### Precision

The precision of the PT determination is highly dependent of the method and technique used. Typical values for within- run and between run precision are:

|          | Within Run (n=10) |        | Between Run (n=10) |        |      |
|----------|-------------------|--------|--------------------|--------|------|
|          | Mean (secs)       | CV (%) | Mean (secs)        | CV (%) |      |
| Sample 1 | 21.3              | 2.36   | Sample 1           | 23.7   | 6.48 |
| Sample 2 | 66.0              | 3.18   | Sample 2           | 67.9   | 3.07 |

**Method Comparison:** Comparison of INR values determined using Helena Capillary Reagent and three other commercially available reagents yielded the following correlation equations:

**Venous Blood**  
Capillary Reagent = 1.1632 (A) - 0.1828 r<sup>2</sup> = 0.9559 n = 20

Capillary Reagent = 0.8317 (B) + 0.1541 r<sup>2</sup> = 0.9593 n = 20

**Capillary Blood**  
Capillary Reagent = 0.9985 (C) - 0.1608 r<sup>2</sup> = 0.9053 n = 29

## BIBLIOGRAPHY

- Quick, A.J. 'A Study of the Coagulation Defect in Hemophilia and Jaundice', Am. J. Med. Sci., 1935; 190: 501.
- Biggs, Rosemary (Ed.), Human Blood Coagulation, Haemostasis and Thrombosis, 2nd Edition, Blackwell Scientific Publications, London, 1976
- Hirsh, J., Poller, L., Deykin, D., Levine, J. and Dalen, J.E., 'Optimal Therapeutic Range for Oral Anticoagulants', Chest, 1989, 95: 55-11S
- Poller, L. 'Laboratory Control of Anticoagulant Therapy', Sem. Thromb. Haemostasis, 1986; 12: 13-19.
- World Health Organisation. Expert Committee on Biological Standards 1984. Technical Series 700:19
- World Health Organisation. Use of Anticoagulants in Diagnostic Laboratory Investigations. 2002; 99.1, Rev.2.
- Poller, L., Triplett, D.A., Hirsh, J., Carroll, J. and Clarke, K. 'The value of plasma calibrants in correcting coagulometer effects on International Normalised Ratios (INR): An international multicentre study' Amer. J. Clin. Pathol., 1995,103 : 358-365.
- Poller, L., Triplett, D.A., Hirsh, J., Carroll, J. and Clarke, K. 'A comparison of lyophilised artificially depleted plasmas and lyophilised plasmas from warfarin treated patients in correcting for coagulometer effects on International Normalised Ratios' Amer. J. Clin. Pathol., 1995, 103 : 366-371.
- The British Society for Haematology, 'Guidelines on Oral Anticoagulation: Second Edition', J. Clin. Pathol., 1990, 43 : 177-183.

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## BIBLIOGRAFIA

- Quick, A.J. 'A Study of the Coagulation Defect in Hemophilia and Jaundice', Am. J. Med. Sci., 1935; 190: 501.
- Biggs, Rosemary (Ed.), Human Blood Coagulation, Haemostasis and Thrombosis, 2nd Edition, Blackwell Scientific Publications, London, 1976
- Hirsh, J., Poller, L., Deykin, D., Levine, J. and Dalen, J.E.,'Optimal Therapeutic Range for Oral Anticoagulants', Chest, 1989, 95: 55-115
- Poller, L. 'Laboratory Control of Anticoagulant Therapy', Sem. Thromb. Haemostasis, 1986; 12: 13-19.
- World Health Organisation. Expert Committee on Biological Standards 1984. Technical Series 700-19.
- World Health Organisation. Use of Anticoagulants in Diagnostic Laboratory Investigations. 2002; 99.1, Rev.2.
- Poller, L., Triplett, D.A., Hirsh, J., Carroll, J. and Clarke, K. 'The value of plasma calibrants in correcting coagulometer effects on International Normalised Ratios (INR): An international multicentre study' Amer. J. Clin. Pathol., 1995; 103: 358-365.
- Poller, L., Triplett, D.A., Hirsh, J., Carroll, J. and Clarke, K. 'A comparison of lyophilised artificially depleted plasmas and lyophilised plasmas from warfarin treated patients in correcting for coagulometer effects on International Normalised Ratios' Amer. J. Clin. Pathol., 1995, 103 : 366-371.
- The British Society for Haematology, 'Guidelines on Oral Anticoagulation: Second Edition', J. Clin. Pathol., 1990; 43 : 177-183.

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- Quick, A.J. 'A Study of the Coagulation Defect in Hemophilia and Jaundice', Am. J. Med. Sci., 1935; 190: 501.
- Biggs, Rosemary (Ed.), Human Blood Coagulation, Haemostasis and Thrombosis, 2nd Edition, Blackwell Scientific Publications, London, 1976
- Hirsh, J., Poller, L., Deykin, D., Levine, J. and Dalen, J.E.,'Optimal Therapeutic Range for Oral Anticoagulants', Chest, 1989, 95: 55-115
- Poller, L. 'Laboratory Control of Anticoagulant Therapy', Sem. Thromb. Haemostasis, 1986; 12: 13-19.
- Poller, L., Triplett, D.A., Hirsh, J., Carroll, J. and Clarke, K. 'The value of plasma calibrants in correcting coagulometer effects on International Normalised Ratios (INR): An international multicentre study' Amer. J. Clin. Pathol., 1995; 103 : 358-365.
- Poller, L., Triplett, D.A., Hirsh, J., Carroll, J. and Clarke, K. 'A comparison of lyophilised artificially depleted plasmas and lyophilised plasmas from warfarin treated patients in correcting for coagulometer effects on International Normalised Ratios' Amer. J. Clin. Pathol., 1995, 103 : 366-371.
- The British Society for Haematology, 'Guidelines on Oral Anticoagulation: Second Edition', J. Clin. Pathol., 1990; 43 : 177-183.

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