Chrom-Z ATIII (anti Xa)

Instructions For Use

REF 5502, 5507

Kit CHROM-Z AT-III (anti-Xa)

Fiche technique

CHROM-Z AT-III (Anti-Xa) Kit

Anleitung

Kit CHROM-Z AT-III (anti-Xa)

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Kit CHROM-Z AT-III (anti Xa)

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Contents

English	1
Français	6
Deutsch	11
Italiano	16
Español	21

INTENDED PURPOSE

The Helena BioSciences CHROM-Z AT-III (anti-Xa) kit is intended for the quantitative determination of antithrombin III (AT-III) activity in human citrated plasma by chromogenic assay. AT-III is a major inhibitor of blood coagulation, acting to inhibit plasma serine proteases including factors IXa, Xa, XIa, and thrombin. The rate of the inhibition is significantly increased in the presence of heparin. AT-III deficiency may be congenital or acquired and is associated with an increased risk for thrombosis¹⁻³. Deficiencies may occur in liver disease, disseminated intravascular coagulation (DIC), septicemia, pulmonary embolism, nephrotic syndrome, stroke, and thrombophlebitis⁴.

Most functional assays for AT-III rely on the ability of plasma to inactivate thrombin in the presence of heparin, such as the assay described by Odegard *et al.*⁴. More recently, it has been shown that a method based on the ability of plasma to inhibit factor Xa in the presence of heparin may be more accurate, since it eliminates interference due to heparin cofactor II⁵ which does not inhibit factor Xa⁶.

In the present two-stage method, factor Xa is added to a plasma dilution containing AT-III in the presence of excess heparin and Calcium. After an initial incubation period (stage 1), residual factor Xa is determined with a factor Xa-specific chromogenic substrate (stage 2). The residual factor Xa activity is inversely proportional to the AT-III concentration.

WARNINGS AND PRECAUTIONS

The reagents contained in this kit are for in-vitro diagnostic use only - **DO NOT INGEST**. Wear gloves when handling all kit components. Refer to the product safety data sheets for risk and safety phrases and disposal information. Plasma products have been screened and found negative (unless otherwise stated on the kit box or vial) for the presence of Hepatitis B Antigen (HbsAg) HIV 1 and 2 antibody and HCV antibody, however they should be handled with the same precautions as a human plasma sample.

Some components of this kit contain sodium azide as a preservative. To prevent the build-up of potentially explosive metallic azides in metal plumbing, flush sinks with copious amounts of water and decontaminate regularly with sodium hydroxide solution.

COMPOSITION**1. Factor Xa Reagent**

Contains freeze-dried bovine factor Xa. Reconstitute each vial with 3x10ml (5502) or 5x2ml (5507) of distilled / deionised water. Allow to stand for 10 minutes and mix well before use.

2. Factor Xa Substrate

Contains freeze-dried CH₃OCO-D-CHA-Gly-Arg-pNA•AcOH. Reconstitute each vial with 3x10ml (5502) or 5x2ml (5507) of distilled / deionised water. Allow to stand for 10 minutes and mix well before use.

3. Sample Diluent:

Contains 4x10ml (5502) or 5x3ml (5507) 5x buffer concentrate with sodium azide at 0.1% as preservative. When fully diluted, buffer contains 0.05 M Tris-HCl, 0.175 M NaCl, 7.5 mM Na₂EDTA and sodium heparin at pH 8.4. Dilute 1 +4 with distilled/deionised water and mix well.

4. Other Kit Components

Each kit contains instructions for use.

STORAGE AND SHELF-LIFE

1. Factor Xa Reagent

Unopened vials are stable until the expiry date indicated on the label. Reconstituted reagent is stable for 2 months at 2...6°C, 1 month at 17°C or 2 days at 37°C. AC-4 POS 32 @ 15°C = 48hrs

2. Factor Xa Substrate

Unopened vials are stable until the expiry date indicated on the label. Reconstituted reagent is stable for 2 months at 2...6°C, 1 month at 17°C or 7 days at 37°C. AC-4 POS 28 @ 37°C = 48hrs.

3. Sample Diluent:

Unopened vials are stable until the expiry date indicated on the label. Store diluted buffer in a tightly stoppered bottle at 2...6°C and use within one month.

ITEMS REQUIRED BUT NOT PROVIDED

(Depending on Assay Method Used):

Cat. No. 5185 SARP

Glacial acetic acid

Coagulation instrument or spectrophotometer operable at 405nm

37°C water bath or dry bath

Laboratory timer

SPECIMEN COLLECTION AND PREPARATION

Plastic or siliconised glass should be used throughout. Blood (9 parts) should be collected into 3.2% or 3.8% sodium citrate anticoagulant (1 part). Separate plasma after centrifugation at 2000xg - 3000xg for 15 minutes. Plasma should be kept at 2...6°C. Testing should be completed within 2 hours of sample collection, or plasma can be stored frozen at -20°C for 2 weeks or -70°C for one month. Thaw quickly at 37°C prior to testing. DO NOT keep at 37°C for more than 5 minutes.

STEP BY STEP PROCEDURE

Preparation of Standard, Control and Patient Dilutions

%AT-III	Plasma	Dilution Buffer
100%	10µl Standard	990µl
50%	500µl 100% Std	500µl
25%	500µl 50% Std	500µl
12.5%	500µl 25% Std	500µl
Patient or Control	10µl Plasma	990µl

Automated Methods, AC-4. (5502 x 400 samples, 5507 x 133 samples)

Refer to the appropriate AC-4 Operators Manual for detailed instructions.

Prepare all reagents as instructed under 'composition'

START=9s	METHOD=Chrom	MIX=0	PATIENT (Pat) vol=4µL pos=35 (CP)	Def. Plasma (DP) Vol=0µL
INCUB=30s	MATH=lin	CLEAN=Y	BUFFER (BUF) Vol=356µL Pos=40	R1(Factor Xa) Vol=75µL Pos=32
RUNTIME=90s	SENS=0	T-CORR=0 S-CORR=0	CLEAR (CLR) Vol=285µL	R2(Substrate) Vol=75µL Pos=28

Semi-Micro Endpoint Method

1. Add 200µl standard / control or patient plasma dilution to a test tube.
2. Incubate at 37°C for 2-4 minutes
3. Add 200µl Factor Xa Reagent and mix.
4. Incubate at 37°C for EXACTLY 1 minute.
5. Add 200µl Factor Xa Substrate and mix.
6. Incubate at 37°C for EXACTLY 3 minutes.
7. Add 200µl acetic acid (50%) and mix.
8. * Add 200µl water (optional)

The yellow colour of the final reaction product is stable for at least 4 hours. Read absorbance at 405nm in a 1cm semi-micro cuvette against a blank prepared in the following order:

1. 200µl acetic acid,
2. 200µl standard dilution,
3. 200µl Factor Xa Reagent,
4. 200µl Factor Xa Substrate,
5. * 200µl water (optional).

* Some spectrophotometers require a minimum of 1ml volume in the cuvette.

If the patient plasma is very icteric, a second blank containing patient plasma dilution instead of standard dilution should be prepared and the absorbance subtracted from the absorbance obtained for the patient's AT-III determination. As many as ten determinations can be performed simultaneously with the same stop watch by staggering pipetting steps at five second intervals.

Kinetic Method

A kinetic analyser may be used to measure the initial rate of hydrolysis of the chromogenic substrate. The procedure to be used is as follows:

To a reaction cuvette:

1. Add 200µl diluted standard or patient plasma.
2. Incubate at 37°C for 1 minutes.
3. Add 200µl Factor Xa Reagent and mix.
4. Incubate at 37°C for 1 minute
5. Add 200µl Factor Xa Substrate
6. Measure rate of change of absorbance at 405nm for 1 minute

Calibration**a) Assay Calibration**

A commercially prepared plasma standard in which AT-III has been determined (eg. Helena BioSciences SARP (Cat. No. 5185)) should be used as a reference.

b) Calibration Curve

Plot the absorbance obtained with each of the AT-III calibration standards on the y-axis against % AT-III on the x-axis using linear graph paper. The line of best fit should be determined by linear regression analysis. The AT-III in plasma samples can be determined by interpolation from the calibration curve. The AT-III concentration in the patient specimen should be adjusted for the AT-III concentration in the standard -

$$\% \text{ AT-III (adjusted)} = \% \text{ AT-III (patient)} \times \% \text{ AT III (standard)} / 100$$

QUALITY CONTROL

Each laboratory should establish a quality control program. Normal and abnormal control plasmas should be tested prior to each batch of patient samples, to ensure satisfactory instrument and operator performance. If controls do not perform as expected, patient results should be considered invalid.

Helena BioSciences supply the following controls available for use with this product:

Cat. No. 5301 SAC-1

Cat. No. 5302 SAC-2

LIMITATIONS

Avoid icteric, lipaemic, and haemolysed samples. To ensure reproducible results, use accurate pipetting devices and observe recommended procedures with emphasis on incubation times and temperature.

REFERENCE VALUES

Reference values can vary between laboratories depending on the techniques and systems in use. For this reason each laboratory should establish its own normal range. The normal range of AT-III is 75-125% in plasma. Plasma activity levels of 30-60% can be observed in patients with hereditary AT-III deficiency^{1,2}. Several clinical conditions associated with acquired AT-III deficiency include: liver disease, DIC, nephrotic syndrome, pulmonary embolism, stroke, and thrombophlebitis. In addition, oral contraceptive use may lead to reduced AT-III levels².

PERFORMANCE CHARACTERISTICS

Helena BioSciences or their representatives have determined the following performance characteristics as a guideline. Each laboratory should establish its own performance data.

Precision

The following estimates of precision (coefficient of variation or C.V.) were observed using an automated analyser utilising a kinetic measurement mode:

Within-assay reproducibility			Between-assay reproducibility		
	AT III Level	CV (%)		AT III Level	CV (%)
n=20	Normal	1.36	n=10	Normal	2.34
n=19	Abnormal	4.40	n=10	Abnormal	4.99

Sensitivity

This assay is designed to give a linear standard curve for AT-III levels between 7.5% and 150%. Samples over 150% should be diluted in Dilution Buffer and re-assayed.

Specificity

The specificity of the assay system has been established in studies employing plasma that has been selectively depleted of AT-III followed by the addition of purified AT-III to achieve various AT-III concentrations.

Accuracy

Comparison of the Helena BioSciences CHROM-Z AT-III kit with another commercially available kit for AT-III determination on 62 samples ranging from 30-130% AT-III activity, the following results were obtained:

$$y = 0.9665x + 0.0629, r^2 = 0.8893$$

BIBLIOGRAPHY

1. Hirsh, J. et al: 'Congenital antithrombin III deficiency. Incidence and clinical features.' *Am J. Med.*, 1989 ; 87 (suppl. 3B) 34S-38S.
2. Panicucci, F. et al: 'Antithrombin III, heparin cofactor and antifactor Xa in relation to age, sex and pathological condition.' *Haemostasis*, 1981; 9: 297-302.
3. Conlan, M.G. et al: 'Antithrombin: Antithrombin III: Associations with age, race, sex and cardiovascular disease risk factors.' *Thromb. Haemost.*, 1994; 72(4) 551-556.
4. Odegard, O.R. et al: 'Heparin cofactor activity measured with an amidolytic method.' *Thromb. Res.*, 1975; 6: 287-294.
5. Demers, C. et al: 'An antithrombin III assay based on factor Xa inhibition provides a more reliable test to identify congenital antithrombin III deficiency than an assay based on thrombin inhibition.' *Thromb. Haemost.*, 1993; 69: 231-235.
6. Tolefsen, D.M., Blank, M.K.: 'Detection of a new heparin-dependent inhibitor of thrombin in human plasma.' *J. Clin. Invest.*, 1981; 68: 589-596.
7. National Committee for the National Laboratory (NCCLS) Standards: Collection transport and preparation of blood specimens for coagulation testing and performance of coagulation assays. Document H21-A2, vol. 11, No. 23, 1991.

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$$y = 0,9665x + 0,0629, r^2 = 0,8893$$

BIBLIOGRAPHIE

1. Hirsh, J. et al: 'Congenital antithrombin III deficiency. Incidence and clinical features'. *Am J. Med.*, 1989 ; 87 (suppl. 3B) 34S-38S.
2. Panicucci, F. et al: 'Antithrombin III, heparin cofactor and antifactor Xa in relation to age, sex and pathological condition.' *Haemostasis*, 1981 ; 9 : 297-302.
3. Conlan, M. G. et al: 'Antithrombin: Antithrombin III: Associations with age, race, sex and cardiovascular disease risk factors.' *Thromb. Haemost.*, 1994 ; 72(4) 551-556.
4. Odegard, O. R. et al: 'Heparin cofactor activity measured with an amidolytic method.' *Thromb. Res.*, 1975 ; 6 : 287-294.
5. Demers, C. et al: 'An antithrombin III assay based on factor Xa inhibition provides a more reliable test to identify congenital antithrombin III deficiency than an assay based on thrombin inhibition.' *Thromb. Haemost.*, 1993 ; 69 : 231-235.
6. Tolefsen, D. M., Blank, M. K.: 'Detection of a new heparin-dependent inhibitor of thrombin in human plasma.' *J. Clin. Invest.*, 1981 ; 68 : 589-596
7. National Committee for the National Laboratory (NCCLS) Standards: Collection transport and preparation of blood specimens for coagulation testing and performance of coagulation assays. Document H21-A2, vol. 11, n° 23, 1991.

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$$y = 0,9665x + 0,0629, r^2 = 0,8893$$

LITERATUR

1. Hirsh, J. et al: 'Congenital antithrombin III deficiency. Incidence and clinical features.' Am J. Med., 1989 ; 87 (suppl. 3B) 34S-38S.
2. Panicucci, F. et al: 'Antithrombin III, heparin cofactor and antifactor Xa in relation to age, sex and pathological condition.' Haemostasis, 1981; 9: 297-302.
3. Conlan, M.G. et al: 'Antithrombin: Antithrombin III: Associations with age, race, sex and cardiovascular disease risk factors.' Thromb. Haemost., 1994; 72(4) 551-556.
4. Odegard, O.R. et al: 'Heparin cofactor activity measured with an amidolytic method.' Thromb. Res., 1975; 6: 287-294.
5. Demers, C. et al: 'An antithrombin III assay based on factor Xa inhibition provides a more reliable test to identify congenital antithrombin III deficiency than an assay based on thrombin inhibition.' Thromb. Haemost., 1993; 69: 231-235.
6. Tolefsen, D.M., Blank, M.K.: 'Detection of a new heparin-dependent inhibitor of thrombin in human plasma.' J. Clin. Invest., 1981; 68: 589-596
7. National Committee for the National Laboratory (NCCLS) Standards: Collection transport and preparation of blood specimens for coagulation testing and performance of coagulation assays. Document H21-A2, vol. 11, No. 23, 1991.

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$$y = 0,9665x + 0,0629, r^2 = 0,8893$$

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1. Hirsh, J. et al: 'Congenital antithrombin III deficiency. Incidence and clinical features.' *Am J. Med.*, 1989; 87 (suppl. 3B) 34S-38S.
2. Panicucci, F. et al: 'Antithrombin III, heparin cofactor and antifactor Xa in relation to age, sex and pathological condition.' *Haemostasis*, 1981; 9: 297-302.
3. Conlan, M.G. et al: 'Antithrombin: Antithrombin III: Associations with age, race, sex and cardiovascular disease risk factors.' *Thromb. Haemost.*, 1994; 72(4) 551-556.
4. Odegard, O.R. et al: 'Heparin cofactor activity measured with an amidolytic method.' *Thromb. Res.*, 1975; 6: 287-294.
5. Demers, C. et al: 'An antithrombin III assay based on factor Xa inhibition provides a more reliable test to identify congenital antithrombin III deficiency than an assay based on thrombin inhibition.' *Thromb. Haemost.*, 1993; 69: 231-235.
6. Tolefsen, D.M., Blank, M.K.: 'Detection of a new heparin-dependent inhibitor of thrombin in human plasma.' *J. Clin. Invest.*, 1981; 68: 589-596
7. National Committee for the National Laboratory (NCCLS) Standards: Collection transport and preparation of blood specimens for coagulation testing and performance of coagulation assays. Document H21-A2, vol. 11, No. 23, 1991.

Helena Biosciences Europe
Queensway South
Team Valley Trading Estate
Gateshead
Tyne and Wear
NE11 0SD

Tel: +44 (0) 191 482 8440
Fax: +44 (0) 191 482 8442
Email: info@helena-biosciences.com
www.helena-biosciences.com

HL-2-1074P 2008/09 (9)