## **Chromogenic Protein C**

Instructions For Use

REF 5543

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**Fiche technique**

**Chromogenes Protein C**

**Anleitung**

**Proteina C Cromogenica**

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**Proteína C Cromógena**

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### **Contents**

|                |    |
|----------------|----|
| English .....  | 1  |
| Français ..... | 6  |
| Deutsch .....  | 12 |
| Italiano ..... | 18 |
| Español .....  | 24 |

**INTENDED USE**

The Helena BioSciences Chromogenic Protein C assay is intended for the quantitative determination of protein C in human plasma using a chromogenic assay method.

Protein C is a vitamin K dependent protein which plays an important role in the regulation of anticoagulant mechanisms. It can inhibit coagulation by inactivating factors Va<sub>1</sub> and VIIIa<sub>2</sub> or, when activated, can stimulate fibrinolysis.<sup>3</sup> Protein C circulates as a zymogen, and is converted to an active serine protease by the action of thrombin in the presence of thrombomodulin. Both hereditary and acquired Protein C deficiencies have been shown to be a risk factor for development of venous thrombosis.

Protein C in plasma is activated by a specific fraction from the Agkistrodon contortrix snake venom. The amount of activated protein C (APC) is determined by monitoring the rate of hydrolysis of a protein C specific chromogenic substrate. The release of pNA is measured at 405 nm and is proportional to the protein C level.

**WARNINGS AND PRECAUTIONS**

The reagents contained in this kit are for In-Vitro diagnostic use only - **DO NOT INGEST**. Wear gloves when handling all kit components. Refer to the product safety data sheets for risk and safety phrases and disposal information.

**COMPOSITION****1. Protein C Substrate**

Each vial contains 2.75  $\mu$ mol lyophilised pyro-Glu-Pro-Arg-pNA.HCl.

**Preparation for Use:** Reconstitute each vial with 2.0ml of Protein C Diluent. If cloudy, warm at 37°C for a few minutes.

**2. Protein C Activator**

Each vial contains 0.8 units of activator from snake venom (Protac®).

**Preparation for Use:** Reconstitute each vial with 2ml deionised water.

**3. Protein C Diluent**

Each vial contains tris buffer with sodium azide as a preservative. Ready for use.

**4. Other kit components**

Each kit contains instructions for use.

**STORAGE AND SHELF LIFE**

Unopened reagents are stable until the given expiry date when stored under conditions indicated on the vial or kit label.

**1. Protein C Substrate**

Reconstituted reagent is stable for 1 week at 2...6°C or one month at -20°C.

AC-4 POS 27@37°C=6hrs

**2. Protein C Activator**

Reconstituted reagent is stable 1 week at 2...6°C or one month at -20°C. AC-4 POS 30@RT = 6hrs.

**3. Protein C Diluent**

Store at 2...6°C.

**ITEMS REQUIRED BUT NOT PROVIDED**

Cat. No. 5185 SARP reference plasma

Coagulation analyser or spectrophotometer capable of reading to 405 nm.

Glacial acetic acid and timer for manual technique

## SAMPLE COLLECTION AND PREPARATION

Plastic or siliconised glass should be used throughout. Blood (9 parts) should be collected into 3.2% or 3.8% sodium citrate anticoagulant (1 part). Separate plasma after centrifugation at 2000-3000 x g for 15 minutes. Plasma should be kept at 2...6°C. Testing should be completed within 2 hours of sample collection, or plasma can be stored frozen at -20°C for 2 weeks or -70°C for one month. Thaw quickly at 37°C prior to testing. Do not keep at 37°C for more than 5 minutes.

Erroneous results may be caused by contamination with tissue fluids or stasis. Avoid agitation, air bubbles or foaming. For the effects of commonly administered drugs, refer to Young, et al. Use only 0.85% NaCl solutions for plasma dilutions.

## STEP BY STEP PROCEDURE

Prepare plasma standards and patient's plasma samples as follows:

| Standard % | Plasma          | Buffer              |
|------------|-----------------|---------------------|
| 100%       | 100 µl S.A.R.P. | + 300µl saline      |
| 50%        | 50 µl S.A.R.P.  | + 350µl saline      |
| 0%         |                 | + 400µl saline only |
| Patient    | 100 µl plasma   | + 300µl saline      |

**IMPORTANT:** Use only 0.85% NaCl solutions for dilutions.

### Automated Methods, AC-4. (5543 x 250 samples)

Refer to the appropriate AC-4 Operators Manual for detailed instructions.

Prepare all reagents as instructed under 'composition'

|              |              |                      |  |   |
|--------------|--------------|----------------------|--|---|
| START=33s    | METHOD=Chrom | MIX=0                | PATIENT (Pat)<br>vol=12µL<br>pos=35 (CP) | Def. Plasma (DP)<br>Vol=0µL             |
| INCUB=300s   | MATH=lin     | CLEAN=Y              | BUFFER (saline)<br>Vol=36µL<br>Pos=39    | PC Activator (R1)<br>Vol=48µL<br>Pos=30 |
| RUNTIME=180s | SENS=0       | T-CORR=0<br>S-CORR=0 | CLEAR (CLR)<br>Vol=0µL                   | PC Substrate (R2)<br>Vol=48µL<br>Pos=27 |

### End Point Method Assay

- To a glass or plastic test tube:  
Add 100 µl standard or patient's plasma dilution  
Incubate at 37°C for 2 minutes  
Add 200 µl Activator and mix.  
Incubate at 37°C for 5 minutes  
Add 200 µl Protein C Substrate and mix  
Incubate at 37°C for 10 minutes  
Add 200 µl acetic acid and mix  
Add (optional) 200µl water
- Read absorbance at 405 nm in a 1cm semi-micro cuvette against a blank prepared with deionised water. Some spectrophotometers require minimum 1ml volume in cuvette. As many as ten determinations can be performed simultaneously with the same stop watch by staggering pipetting steps at five second intervals.

**Kinetic Method**

Automated analyser required. Refer to Operators Manual for detailed instructions.

Example table:


| Reagent             | Volume     | Time |
|---------------------|------------|------|
| Sample              | 75 $\mu$ l | 3:00 |
| Activator           | 75 $\mu$ l | 5:00 |
| Protein C Substrate | 75 $\mu$ l | 5:00 |

**Calibration Curve**

Plot the absorbance obtained for each of the protein C standards against protein C % on linear graph paper.

The protein C concentration in patient's plasma specimens can be determined by interpolation from the calibration curve. If a commercial protein C standard is used, the protein C concentration in the patient's specimen should be adjusted for the protein C concentration in the standard.

The calibration curve shown below is an example only. A calibration curve must be obtained each time the assay is performed.

**INTERPRETATION OF RESULTS**

Protein C deficiencies, either congenital or acquired, may lead to serious thrombotic events such as thrombophlebitis, deep vein thrombosis, or pulmonary embolism. Approximately 2-8% of all patients with venous thrombosis under the age of 40-45 years have Protein C deficiencies. Patients with hereditary deficiencies generally present with venous thrombosis in young adulthood. The first episode is usually spontaneous and associated with trauma or stress to the haemostasis mechanism. The prevalence of Protein C deficiencies is one case per 300 people or approximately 0.33%.

## Acquired Deficiencies

Decreased levels of protein C are observed in the following cases:

- hepatic disorders: hepatitis, cirrhosis.
- DIC
- oral anticoagulant therapies, in which cases the interpretation of test results is difficult if the patients have had a history of thromboses and are receiving oral anticoagulant treatments.

## Congenital Deficiencies

Two types of deficiencies can be observed:

- quantitative type I, in which the functional and antigenic levels are simultaneously decreased
- qualitative or type II, characterised by a decreased functional level and by a normal or subnormal immunological level.

## QUALITY CONTROL

Each laboratory should establish a quality control program. Normal and abnormal control plasmas should be tested prior to each batch of patient samples, to ensure satisfactory instrument and operator performance. If controls do not perform as expected, patient results should be considered invalid.

Helena BioSciences supply the following controls available for use with this product:

Cat. No. 5301 SAC 1

Cat. No. 5302 SAC 2

## LIMITATIONS

Icteric, haemolysed or hyperlipemic specimens interfere with absorbance readings, thus requiring the use of plasma blanks for accurate results. Plasma blanks are also needed on patients when contact factor activation is suspected, such as DIC patients, or from individuals on oral contraceptives where cold activation may occur. Blanks are done by substituting saline for the Activator in the test reaction. Subtract the sample blank activity from the sample test activity.

The naturally-occurring activator of protein C is thrombin in the presence of thrombomodulin. The possible existence of clinical protein C abnormalities not detectable by reactivity with the snake venom activator or with the chromogenic substrate used in this kit should not be precluded. To assure accurate, reproducible results, use accurate pipetting devices and observe recommended procedures with particular emphasis on incubation times and incubation temperature.

## REFERENCE VALUES

Reference values can vary between laboratories depending on the techniques and systems in use. For this reason each laboratory should establish its own normal range. Protein C activity values are usually expressed in relative percentages compared to a pooled normal plasma standard. Bertina, et al.<sup>5</sup> reported a range of 65-145% in healthy individuals with diminished levels found following anticoagulant therapy. There is apparently no difference in protein C between healthy males and females.<sup>6</sup> Helena BioSciences tested 17 plasmas of presumed healthy men and women using the Helena AC-4 with the following data.

|                     |
|---------------------|
| X = 105%            |
| SD = 26.6           |
| range = 61.2-162.9% |

**PERFORMANCE CHARACTERISTICS**

Helena BioSciences or their representatives have determined the following performance characteristics as a guideline. Each laboratory should establish its own performance data.

**I. Specificity**

In studies where recovery of protein C plasma was determined after addition of various amounts of protein C to protein C deficient plasma, the expected recovery was obtained. **II.**

**Precision**

Within Run: Precision was run on a control in replicate, resulting in the following data:

|           |
|-----------|
| n = 10    |
| X = 68.5  |
| SD = 1.7  |
| CV% = 2.5 |

Between Run: Precision was run on a control in replicate. The data from 3 consecutive runs was combined, and gave the following results:

|           |
|-----------|
| n = 18    |
| X = 45.3  |
| SD = 1.2  |
| CV% = 2.6 |

**BIBLIOGRAPHY**

1. Walker, F.J., J Biol Chem 256:11128-11131, 1981.
2. Taylor, F.B. and Lockhart, M.S., Thrombin Res 27:115-164, 1985.
3. Fulcher, C.A., et al., Blood 63:486-489, 1981.
4. Young, D.S. et al., Effects of Drugs on Clinical Laboratory Tests, 3rd ed., AACC Press, Washington, D.C., 1990
5. Bertina, R.M., Broekmans, A.W., Krommenhoek-Van Es, C., Van Wijngaarden, A. The Use of a Functional and Immunologic Assay for Plasma Protein C in the Study of the Heterogeneity of Congenital Protein C Deficiency. Thromb Haemostas, 51:1-5, 1984.
6. Pabinger-Fasching, I., Bertina, R.M., Lechner, K., Niesser, H., Korninger, Ch. Protein C Deficiency in Two Austrian Families. Throm Haemostas, 50:810-813, 1983.

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|  |
|--|
| $n = 18$<br>$X = 45,3$<br>$SD = 1,2$<br>$CV\% = 2,6$ |
|--|

**BIBLIOGRAFIA**

1. Walker, F.J., J Biol Chem 256:11128-11131, 1981.
2. Taylor, F.B. and Lockhart, M.S., Thrombin Res 27:115-164, 1985.
3. Fulcher, C.A., et al., Blood 63:486-489.
4. Young, D.S. et al., Effects of Drugs on Clinical Laboratory Tests, 3rd ed., AACC Press, Washington D.C., 1990
5. Bertina, R.M., Broekmans, A.W., Krommenhoek-Van Es, C., Van Wijngaarden, A. The Use of a Functional and Immunologic Assay for Plasma Protein C in the Study of the Heterogeneity of Congenital Protein C Deficiency. Thromb Haemostas, 51:1-5, 1984.
6. Pabinger-Fasching, I., Bertina, R.M., Lechner, K., Niesser, H., Korninger, Ch. Protein C Deficiency in Two Austrian Families. Throm Haemostas, 50:810-813, 1983.

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