INTENDED PURPOSE

The Fibrinogen Assay Kit is intended for the quantitative determination of fibrinogen in human plasma.

Clauss' developed a simple method for the quantitative determination of fibrinogen by measuring the clotting time of dilute plasma after the addition of thrombin (>30 NIH units/ml). This clot time is proportional to the fibrinogen concentration. Levels of fibrinogen can increase as a result of inflammation, pregnancy or oral contraceptive use⁶. Decreased levels can be found in certain states such as liver disease and DIC. Congenital deficiencies include afibrinogenaemia (no detectable fibrinogen), hypofibrinogenaemia (<1 mg/ml) and dysfibrinogenaemia (abnormal fibrinogen molecule).

WARNINGS AND PRECAUTIONS

The reagents contained in this kit are for in-vitro diagnostic use only. **DO NOT INGEST.** Wear gloves when handling all kit components. Refer to the product safety data sheets for risk and safety phrases and disposal information. Plasma products have been screened and found negative (unless otherwise stated on the kit box or vial) for the presence of Hepatitis B Antigen (HbsAg) HIV 1 and 2 antibody and HCV antibody; however they should be handled with the same precautions as a human plasma sample.

COMPOSITION

- Thrombin Reagent (5 x 4ml)**
Contains approximately 200 units bovine thrombin with stabilisers.
Preparation: Reconstitute each vial with 4.0ml of purified water. Take care when pipetting to avoid contamination. Swirl gently and allow to stand for 15 minutes. Mix well before use. Do not shake.
- Fibrinogen Calibrator (2 x 1ml)**
Contains 1ml of lyophilised normal human plasma assayed for fibrinogen levels.
Preparation: Reconstitute each vial with exactly 1.0ml of purified water. Swirl gently and allow to stand for 10 minutes. Mix gently before use. Do not shake.
- Owren's Veronal Buffer (2 x 25ml)**
Contains 25ml of buffer which contains barbital and sodium chloride and sodium azide as preservative. The buffer is ready for use as packaged.
- Other kit components**
Each kit contains Instructions For Use and lot specific reference values insert.

STORAGE AND SHELF-LIFE

Unopened vials are stable until the given expiry date when stored under conditions indicated on the vial or kit label.

- Thrombin Reagent**
Once reconstituted, the reagent is stable for 8 hours at 15...30°C or 1 week at 2...6°C, AC-4 POS 29 @RT = 8hrs.
- Fibrinogen Calibrator**
Once reconstituted, the reagent is stable for 4 hours at 2...6°C.
- Owren's Veronal Buffer**
Store at 2...6°C once opened, AC-4 POS 39 @37°C = 8hrs.

SAMPLE COLLECTION AND PREPARATION

Plastic or siliconised glass should be used throughout. Blood (9 parts) should be collected into 3,2% or 3,8% sodium citrate anticoagulant (1 part). Separate plasma after centrifugation at 2000 - 3000 x g for 15 minutes. Plasma should be kept at 2...6°C. Testing should be completed within 2 hours of sample collection, or plasma can be stored frozen at -20°C for 2 weeks or -70°C for one month. Thaw quickly at 37°C prior to testing. Do not keep at 37°C for more than 5 minutes.

STEP BY STEP PROCEDURE

Helena BioSciences Fibrinogen Assay reagents are suitable for use on any coagulation analyser capable of performing Clauss Fibrinogen assays. Refer to the Instrument Operators Manual and/or application note for appropriate instructions.

- A. Automated Methods, AC-4. (x400 samples)**
Refer to the appropriate AC-4 Operators Manual for detailed instructions. Prepare all reagents as instructed under 'composition'.

START =3s	METHOD =Coag	MIX=0	PATIENT (Pat) vol=10µL Pos=35 (CP)	Def. (DP) Vol=0µL	Plasma
INCUB =120s	MATH =logXY	CLEAN=Y	BUFFER (BUF) Vol=90µL Pos=39	R1 Vol=0µL	
RUNTIME =120s	SENS =0	U-CORR=0 S-CORR=0	CLEAR (CLR) Vol=0µL	R2 Vol=50µL Pos=29	

QUALITY CONTROL

Each laboratory should establish a quality control program. Normal and abnormal control plasmas should be tested prior to each batch of patient samples, to ensure satisfactory instrument and operator performance. If controls do not perform as expected, patient results should be considered invalid.

- Helena BioSciences supply the following controls available for use with this product:
 REF 5301 Sac-1 Normal Control Plasma 10x1ml
 REF 5302 Sac-2 Abnormal Control Plasma 10x1ml
 REF 5186 Norm-Trol 1 Control Plasma 10x1ml
 REF 5187 Ab-Trol 2 Control Plasma 10x1ml
 REF 5183 Ab-Trol 3 Control Plasma 10x1ml

REFERENCE VALUES

Reference values can vary between laboratories depending on the techniques and systems in use. For this reason each laboratory should establish it's own normal range. Expected values for fibrinogen in healthy adults are 150-350 mg/dl (1.5-3.5g/L)3,4.

LIMITATIONS

Heparin levels >0.6 units/ml and fibrinolytic degradation products >100mg/ml may cause falsely low fibrinogen quantitation. If values fall outside the standard curve values for the patient samples, re-assay using an appropriate dilution to bring values into the standard range.

PERFORMANCE CHARACTERISTICS

Helena BioSciences or their representatives have determined the following performance characteristics as a guideline. Each laboratory should establish it's own performance data. The Fibrinogen assay is designed to give a linear calibration from 1.5 - 6.5g/L.

REPRODUCIBILITY:

Within assay reproducibility:		
Fibrinogen (g/L)	CV	n
1.24	4.88	5
3.01	3.58	5

Between assay reproducibility:		
Fibrinogen (g/L)	CV	n
1.02	6.17	10
3.62	3.78	10

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