Instructions For Use

APTT-SiL Minus

REF 5558, 5559

APTT-SiL Minus

Fiche technique

APTT-SiL Minus

Anleitung

APTT-SiL Minus

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APTT-SiL Minus

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Contents

English	1
Français	6
Deutsch	11
Italiano	16
Español	21

INTENDED USE

The Activated Partial Thromboplastin Time (APTT) test is used principally to detect deficiencies in the Stage I Factors, namely Factors VIII, IX, XI and XII. From its origins through the work of Langdell and coworkers¹, and later modified by others^{2,3}, the Activated Partial Thromboplastin Time test has been widely used for a number of years as a pre-surgical screen for assessing certain coagulation factors and in monitoring Heparin therapy. All Factors of the Intrinsic Pathway are necessary for normal results when performing the APTT test.

The APTT test is also used to monitor Heparin (UFH) therapy, showing prolonged test results at approximately 0.1 units/ml and above⁴.

The test is also used in the quantitative determination (Factor Assays) of Factors VIII, IX, XI, XII.

The APTT test measures the clotting time of test plasma. An initial activation of the APTT reagent with the test plasma is followed by the addition of calcium chloride and the clotting time is measured in seconds. Deficiencies of approximately 40% and lower of Factors VIII, IX, XI and XII will result in a prolonged APTT. Heparin, in the presence of adequate amounts of AT-III will also result in a prolonged APTT.

WARNINGS AND PRECAUTIONS

All wastes containing biological material should be properly labelled and stored separately from other wastes. Dispose of all waste materials according to prescribed international, national and local regulations.

The test should be used in conjunction with clinical observations and results of other laboratory tests.

COMPOSITION

REF 5558 (5 x 10ml APTT-SiL Minus Reagent + 5 x 10ml 25mM Calcium Chloride)

REF 5559 (10 x 10ml APTT-SiL Minus Reagent + 10 x 10ml 25mM Calcium Chloride)

APTT-SiL Minus reagent contains a near-collodial particle activator (magnesium-aluminum-silicate) for optimum sensitivity to Factor deficiencies and to Heparin. The Reagent also contains phospholipids with buffer and stabilizers.

Preparation for use:

Bring to room temperature prior to use. Mix well by swirling or inversion prior to use.

STORAGE AND STABILITY

Store at 2...8°C. Do not freeze.

Stable for 30 days after opening. Avoid prolonged heating.

SAMPLE COLLECTION AND PREPARATION

- A. Anticoagulant. Use buffered sodium citrate 3.2%.
- B. Specimen Collection⁶
 1. Obtain venous blood by clean venipuncture.
 2. Immediately mix nine parts blood with one part anticoagulant, mix well by inversion of tube against the stopper.
 3. Centrifuge the specimen at 1500 x g for no less than 15 minutes.
 4. Remove plasma from the tube within 60 min using a plastic pipette and store in a plastic tube.
 5. Test plasma sample within 4 hours, otherwise store frozen at -20°C for up to two weeks or -70°C for up to six months⁶ and thaw once just prior to use.

STEP-BY-STEP PROCEDURE

A. Materials Provided:

1. APTT-SiL Minus (liquid)
2. 0.025 M Calcium chloride

B. Materials Required:

1. A manual, mechanical or photo-optical means of clot detection.
2. Suitable timer.
3. Parafilm or equivalent.
4. Precision pipettor, 0.1 ml.
5. Control Plasmas.

C. Performance of Test

The use of APTT-SiL Minus in the performance of APTT determinations is adaptive to manual or automated instruments for clot detection. See instrument manufacturers instructions for full details.

1. Prewarm well mixed APTT-SiL Minus and 0.025 M Calcium Chloride to 37°C.
2. Prewarm 0.1 ml test plasma in duplicate to 37°C for 2 minutes.
3. Forcibly add 0.1 ml prewarmed APTT-SiL Minus to plasma and start timer.
4. Incubate for exactly 5 minutes at 37°C.
5. Add 0.1 ml prewarmed 0.025 M Calcium Chloride.
6. Note time for clot formation.
7. Report result as APTT time (seconds).

D. Automated Methods, AC-4. (5558 x1000 samples, 5559 x2000 samples)

Refer to the appropriate AC-4 Operators Manual for detailed instructions.

Prepare all reagents as instructed under 'composition'.

START=17s	METHOD=Coag	MIX=0	PATIENT(Pat) vol=50µL pos=35 (CP)	Def. Plasma (DP) Vol=0µL
INCUB=300s	MATH=lin	CLEAN=0	BUFFER (BUF) Vol=0µL	R1Vol=50µL Pos=31
RUNTIME=180s	SENS=I	T-CORR=0 S-CORR=0	CLEAR (CLR) Vol=0µL	R2Vol=50µL Pos=26

INTERPRETATION OF RESULTS

Calculate the mean clotting time of duplicate determinations for each test plasma and Control Plasmas. The difference between duplicate determinations should not exceed 5% (repeat if necessary). Report the APTT time result in seconds along with the Normal Range.

Control Plasmas are to be used to help assure that test variables, such as temperature, reagents, instrumentation, pipettes, etc. are functioning properly. It is not recommended that Control values be reported with patient results.

QUALITY CONTROL

Helena Normal and Abnormal Coagulation Control Plasmas should be tested in conjunction with patient samples. It is recommended that at least one Normal and one Abnormal be run at each shift.

A Control Range should be established by the laboratory to determine the allowable variation in day to day performance of each Control Plasma.

EXPECTED VALUES

Typical Normal results using Helena APTT-SiL Minus are approximately 30-40 seconds. However, each laboratory should establish its own Normal Range using the reagents, instrumentation and techniques commonly in use in the laboratory. A Normal Range should be established using individuals representative of the population being tested within the laboratory.

A new Normal Range should be established with any change in instrumentation, blood collection techniques, anticoagulant and when changing to new lots of reagents.

LIMITATIONS

Great care must be taken to minimize variations which may occur by seemingly insignificant factors.

A. Specimen Collection.**AVOID:**

1. Use of glass. Use only plastic.
2. Delayed mixing of blood with anticoagulant.
3. Hemolysis or lipemic samples.
4. Contamination with tissue thromboplastin.
5. Improper ratio of anticoagulant with blood.

B. Laboratory Techniques

1. Perform tests at 37°C.
2. Use only high purity water.
3. Optimum pH is 7.0-7.5

PERFORMANCE CHARACTERISTICS

Helena BioSciences or their representatives have determined the following performance characteristics as a guideline. Each laboratory should establish its own performance data.

Precision Studies

Results are shown below:

Within Run (aPTT in seconds): Precision studies were performed in replicate on the Sysmex CA1500 using two levels of control materials.

Control	Mean	C.V.	N
Norm-Trol 1	33	0.36%	10
Ab-Trol 3	77.9	0.31%	10

Between Run (aPTT in seconds): Precision studies were performed in the same manner using two levels of control material during three runs for evaluation.

Control	Mean	C.V.	N
Norm-Trol 1	32.9	2.41%	100
Ab-Trol 3	78.3	0.77%	100

Factor Sensitivity

Studies were performed to assess the sensitivity of APTT-SiL Minus reagent to various Factor Deficiencies. The following results were obtained.

% Factor	Factor VIII (Seconds)	Factor IX (Seconds)	Factor XI (Seconds)
<1	85.7	70.9	92.4
10	45.9	44.4	50.7
40	34.2	33.9	34.4
100	28.8	28.8	28.8

Upper limit of Normal for this study = 34.9 seconds.

Heparin Sensitivity

Numerous variables can affect the behaviour of heparinized plasma samples when evaluated by the APTT test. As a general guideline, when increasing amounts of heparin were added to normal plasma, the following results were obtained using APTT-SiL Minus on a photo-optical instrument.

Heparin IU/ml	Clot Time (seconds.)
0	29.9
0.2	70.0
0.4	174

REFERENCES

1. Langdell, R., Wagner, R., Brinkhous, K.: Effect of Antihemophilic Factor on One Stage Clotting Tests. *J. Lab. Clin. Med.* 41: 637, 1953.
2. Proctor, R., Rapaport, S.: The Partial Thromboplastin Time with Kaolin. *Am. J. Clin. Path.* 36: 212, 1961.
3. Ratnoff, O., Crum, S.: Activation of Hageman Factor by Solutions of Ellagic Acid. *J. Lab. Clin. Med.* 63: 359, 1964.
4. Hathaway, W., Clarke, S., Christian, M., Jacobson, L.: Sensitivity of the Activated Partial Thromboplastin Time to Mild Factor Deficiencies; in *Standardization of Coagulation Assays: An Overview*. College of American Pathologists, Skokie, IL. pp. 131, 1982.
5. Brandt, J., Triplett, D.: Laboratory Monitoring of Heparin. Effect of Reagents and Instruments in the Activated partial Thromboplastin Time. *Am. J. Clin. Path.* 76: 530, 1981.
6. CLSI: Guidelines for the Standardized Collection, Transport and Preparation of Blood Specimens for Coagulation Testing and Performance of Coagulation Assays. H21-A4.

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LITERATUR

1. Langdell, R., Wagner, R., Brinkhous, K.: Effect of Antihemophilic Factor on One Stage Clotting Tests. J. Lab. Clin. Med. 41: 637, 1953.
2. Proctor, R., Rapaport, S.: The Partial Thromboplastin Time with Kaolin. Am. J. Clin. Path. 36: 212, 1961.
3. Ratnoff, O., Crum, S.: Activation of Hageman Factor by Solutions of Ellagic Acid. J. Lab. Clin. Med. 63: 359, 1964.
4. Hathaway, W., Clarke, S., Christian, M., Jacobson, L.: Sensitivity of the Activated Partial Thromboplastin Time to Mild Factor Deficiencies; in Standardization of Coagulation Assays: An Overview. College of American Pathologists, Skokie, IL. pp. 131, 1982.
5. Brandt, J., Triplett, D.: Laboratory Monitoring of Heparin. Effect of Reagents and Instruments in the Activated partial Thromboplastin Time. Am. J. Clin. Path. 76: 530, 1981.
6. CLSI: Guidelines for the Standardized Collection, Transport and Preparation of Blood Specimens for Coagulation Testing and Performance of Coagulation Assays. H21-A4.

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