Potassium-p

Presentation: 2 x 50 mL. Ref.: 30343 Store: 2 - 8 °C.

POTASSIUM TPB-Na Method

Quantitative determination of Potassium.

Only for *in vitro* use in clinical laboratory (IVD)

TEST SUMMARY

Potassium ions in a protein-free alkaline medium react with sodium tetraphenylboron to produce a finely dispersed turbid suspension of potassium tetraphenylboron. The turbidity produced is proportional to the potassium concentration and read photometrically.

REAGENTS COMPOSITION

PREC	Trichloroacetic acid (TCA)	0.3 mol/L
R.1 TPB-Na	Sodium tetraphenylboron (TPB-Na)	0.2 mol/L
R.2 NaOH	Sodium hydroxide	2.0 mol/L
K-p CAL	Sodium aqueous primary standard 5.0 mmol/L	

PRECAUTIONS

R.2 (naoh): C: R35 Irritating to eyes and skin. PREC, CAL (TCA) : Xi R 36/38 Irritating to eyes and skin. S24/25 avoid contact with skin and eyes. S26 in case of contact with eyes, rinse immediately with plenty of water and seek medical advice. S37/39: wear suitable gloves and eye/face protection. S45: in case of accident or if you feel unwell, seek medical advice immediately (show the label where possible).

REAGENT PREPARATION AND STABILITY

Working reagent (wr): Mix equal volumes of R.1 TPB-Na and R.2 naoh. Allow to stand for 15-30 minutes prior to use. The working reagent is stable for 7 days at 15-25°C and 30 days at 2-8°C.

All the reagents of the kit are stable up to the end of the indicated month and year of expiry. Store tightly closed at 2-8°C. Do not use reagents over the expiration date.

SPECIMEN

Non-haemolytic serum or heparin plasma

MATERIAL REQUIRED BUT NOT PROVIDED

- Spectrophotometer or colorimeter measuring at 578 nm.
- Matched cuvettes 1.0 cm light path.

General laboratory equipment.

TEST PROCEDURE

- Assay conditions:
 - Wavelength: 578 nm
 - Cuvette: 1 cm. light path
 - Temperature: 37°C /15-25°C
- Adjust the instrument to zero with distilled water.
- Pipette into a cuvette:

Sample (µL)	25
Precipitating sol. (µL)	250

- Mix carefully.
- Centrifuge at high speed for 5-10 min.
- Separate the clear supernatant and pipette on another cuvette:

	Standard	Sample
Working reagent (µL)	500	500
Standard (µL)	50	--
Supernatant (µL)	--	50

- To produce an homogeneous turbidity, the standard or the clear supernatant must be added to the center of the surface of the working reagent in the cuvette. Mix each cuvette carefully before proceeding to the next sample. Mix and allow for stand for 5 min.
- Read the absorbance (A) of standard and samples against working reagent blank between 5 and 30 minutes.

CALCULATIONS (Note 2)

$$\frac{A_{\text{Sample}}}{A_{\text{STD}}} \times 5.00 \text{ (Standard conc.)} = \text{mmol/L potassium in the sample}$$

Conversion factor: mmol/L = mEq/L.

QUALITY CONTROL

Control sera are recommended to monitor the performance of the procedure, LABTROL H Normal Ref. 30950 and LABTROL H Pathological Ref. 30955. If control values are found outside the defined range, check the instrument, reagents and calibrator for problems.

Serum controls are recommended for internal quality control. Each laboratory should establish its own Quality Control scheme and corrective actions.

REFERENCE VALUES

Serum:	3.60 – 5.50 mmol/L
Plasma:	4.00 – 4.80 mmol/L

It is suggested that each laboratory establish its own reference range.

CLINICAL SIGNIFICANCE

Potassium (K+) is the major positive ion within cells and is particularly important for maintaining the electric charge on the cell membrane. This charge allows nerves and muscles to communicate and is necessary for transporting nutrients into cells and waste products out of the cell. The concentration of potassium inside cells is about 30 times that in the blood and other fluids outside of cells. Potassium levels are mainly controlled by the steroid hormone aldosterone. Aldosterone is secreted from the adrenal gland when levels of potassium increase. Aldosterone, in turn, causes the body to rid itself of the excess potassium. Metabolic acidosis (for example, caused by uncontrolled diabetes) or alkalosis (for example, caused by excess vomiting) can affect blood potassium. In normal people, taking potassium supplements or potassium-containing drugs is of no consequences, because the kidneys efficiently dispose of excess potassium.

REAGENT PERFORMANCE

1. **Measuring range:** From detection limit of 2 mmol/L to linearity limit of 20 mmol/L.

If the results obtained were greater than linearity limit, dilute the sample 1/2 with NaCl 9 g/L and multiply the result by 2.

2. **Precision:**

	Intra-assay (n=20)		Inter-assay (n=20)	
	Mean (mmol/L)	SD	Mean (mmol/L)	SD
Mean (mmol/L)	4.15	6.70	4.15	6.70
SD	0.11	0.176	0.152	0.19
CV (%)	2.58	2.54	4.11	2.23

3. **Sensitivity:** 1 mmol/L = 0.537A.

4. **Accuracy:** Results obtained using LABKIT reagents did not show systematic differences when compared with other commercial reagents.

The results of the performance characteristics depend on the analyzer used.

INTERFERING SUBSTANCES

A list of drugs and other interfering substances with potassium determination has been reported by Young et. al^{5,6}.

NOTES

- K-p CAL: Proceed carefully with this product because due its nature it can get contaminated easily.
- As red blood cells contain about 25 times the amount of potassium, they have to be separated from the serum within one hour after blood collection. Otherwise, falsely elevated potassium concentrations will be found.
- Traces of detergents produce turbidity which leads to falsely elevated potassium concentrations. They therefore have to be avoided.
- Calibration with the aqueous standard may cause a systematic error in automatic procedures. In these cases, it is recommended to use a serum Calibrator.
- CHEMELEX has instruction sheets for several automatic analyzers. Instructions for many of them are available on request.

BIBLIOGRAPHY

- Hillmann, G., Beyer, G., Z. Klin. Chem. Klin. Biochem. 5, 93 (1967)
- Henry, R.J., Clin. Chem., Harper & Row, New York, Sec. Edit. 646 (1974)
- Tietz, N.W., Fundamentals of Clinical Chemistry, Saunders, Philadelphia, Sec. Edit., 876 (1976)
- ISO 15223 Medical devices – Symbols to be used with medical device labels, labelling and information to be supplied.
- Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
- Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.

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BIBLIOGRAFIA

- 1. Hillmann, G., Beyer, G., Z. Klin. Chem. Klin. Biochem. 5, 93 (1967)
2. Henry, R.J., Clin. Chem., Harper & Row, New York, Sec. Edit. 646 (1974)
3. Tietz, N.W., Fundamentals of Clinical Chemistry, Saunders, Philadelphia, Sec. Edit., 876 (1976)
4. ISO 15223 Medical devices - Symbols to be used with medical device labels, labelling and information to be supplied.
5. Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
6. Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.